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CE 0197

F - 01/2010

HOSEPIPE FOR LOW-PRESSURE MEDICAL FLUIDS

In compliance with Standards NF EN ISO 5359

Checks before use: if the tube is dropped, the mistake-proof notch may be damaged, making connection impossible. Always check your set of tubes before use.


- Safety:**
- Do not use close to any ignition point
 - Protect the tube against exposure to high temperatures
 - Do not bring the tube into contact with greasy products
 - These tubes should only be used by staff members who have received training on the hazards associated with using medical gases.

Device: low-pressure tube with gas-specific nozzles; in compliance with Standard NFS 90-116. The nozzles are permanently and irreversibly fixed onto the tube. Device designed for connection to standard outlets on pipes for compressed gases.

Connection:

- **First**, connect the straight nozzle of the tube to the device to be used and **then** connect the T-shaped nozzle to the tank, pressing on the nozzle and turning it to the right.

Expiration and traceability

- Device lifespan: 10 years
- Shown on a pictogram on the T-shaped nozzle
- E.g.:  (change before 2020)
- Engraved on the T-shaped nozzle: month and year of manufacture together with the Validation Number TM CE 0197
- Tube Batch N° (ink-jet-printed)
- Label on the packaging with the tube expiration date and Production Batch N°

Maintenance and checking:

- An annual inspection (minimum) should be carried out
- Overall condition
 - Airtightness
 - Inspecting the nozzles: check the notches on the nozzle and make sure that it cannot be connected to a tank containing a different gas

Specifications:

- Class II a device
 - Delivered NON-sterile
 - Maximum pressure: 15 bar
 - Bursting pressure: 60 bar
 - Operating temperature: - 15 to + 60°C
- Inner diameter: 6.0 mm \pm 0.3
Outer diameter: 12.00 mm \pm 0.3
Contains no latex
Contains phthalates

Cleaning:

- Clean with a regular detergent in water
- Do not immerse the tube

Storage:


- keep away from ultraviolet radiation at between - 10°C and + 40°C

The TM tube for medical gases is made of anti-static PVC in four layers: extra-flexible with reinforcement and high-resistance polyester fibre.

- Odourless
- Smooth with a glaze finish, very easy to keep clean
- Highly resistant to compression and puncture
- The PVC is flexible and light

In line with European Directive 93/42/CEE of 14 June 93, any incident or risk of an incident must be reported to Medical Technology immediately.

- Technical characteristics:**
- Class IIa device
 - Delivered non-sterile
 - Lifespan: 10 years
 - Expiration date shown on nozzle
 - Operating pressure: up to 15 bar
 - Bursting pressure: 60 bar
 - Operating temperature: - 15 to + 60°C
 - Inner diameter: 6.0 mm \pm 0.3
 - Outer diameter: 12.00 mm \pm 0.3

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