DETENDEURS AVEC DEBITMETRE A ORIFICES PRECALIBRES
PRESSURE REGULATORS WITH FLOWMETER WITH PRE-ADJUSTED FLOWRATES

DETENDEURS REGSON TM2

REGSON TM2 PRESSURE REGULATORS

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Medical device meant to be connected to high pressure (O₂ or medical Air) medical gas cylinders, to regulate the gas pressure by supplying the operating outlet with a controlled and stable pressure below its supply pressure. The gauge indicates the pressure level of the gas cylinder. It is combined with a DEBSON TM2 flowmeter with pre-adjusted flowrates.

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• **Classification rule and CE certificate registration No.:**

- Active device of class IIb non sterile - according to the classification of the annexe IX of the European directive 93/42/EEC and 2007/47/EEC, chapter 3 rule 11.
- CE 0197 TÜV RHEINLAND
- HD 60020766 0001
- Directive 93/42/EEC and 2007/47/EEC, Annexe V, Article 3

• **UMDNS code and product group:**

11748: "Flowmeter gas"

• **GMDN code and product group:**

36271: "A system designed to supply compressed medical gases (oxygen, nitrous oxide, air) and vacuum from a source to outlets and inlets throughout a medical facility. A system may consist of a gas supply, manifold unit, a pressure regulation and alarm system, a pipework system, remote or local indicators, pressure relief valves, zone valves, and wall outlet terminal units."

• **Applied standards to the products:**

EN ISO 10524-1: 2006 - Pressure regulators for use with medical gases. Part 1: pressure regulators and pressure regulators with flow-metering devices.

EN ISO 15002: 2008 - Flow-metering devices for connection to terminal units of medical gas pipeline systems.

DETENDEURS AVEC DEBITMETRE A ORIFICES PRECALIBRES
PRESSURE REGULATORS WITH FLOWMETER WITH PRE-ADJUSTED FLOWRATES

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- **Available gases**: Oxygen (O₂) and Medical Air

- **Supply pressure**: Up to 200 bar maximum

- **Working pressure**: 4.5 bar +/-0.5

- **Unit of measurement**:

- Pressure: bar
- Flowrate: liter per minute (l/min)

- **Reading accuracy**: ± 10 % of the read value

- **Working temperature**: 23°C

- **Flowrates breakdown**:

01 l/min:	0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	1
05 l/min:	0	0.5	0.8	1	1.5	2	2.5	3	4	5
15 l/min:	0	1	2	3	4	5	6	9	12	15
25 l/min:	0	0.5	1	2	3	4	6	8	15	25
50 l/min:	0	5	10	15	20	25	30	35	40	50

- **Connection systems to the gas cylinder**:

- **Bullnose fitting**, for large cylinder (> 4 L in water or 0.8 m³) according to various standards: AFNOR (French Standard), BS (British Standard), DIN (German Standard), US OHMEDA DIAMOND (American Standard), NORDIC (Scandinavian Standard), UNI (Italian Standard) etc.
- **Yoke pin-index fitting (universal system)**, for small cylinder.

- **Available models**:



- Available with or without quick-release connector
- Available with or without Venturi ejector system (OXAP)
- Available with or without Flush system (for 15 l/min only)

- **Weight**:

- With AFNOR bullnose fitting: 828 g.
- With universal pin-index fitting: 755 g.

- **Dimensions** (with AFNOR bullnose fitting):

DETENDEURS AVEC DEBITMETRE A ORIFICES PRECALIBRES
PRESSURE REGULATORS WITH FLOWMETER WITH PRE-ADJUSTED FLOWRATES

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• **Technical data:**

Piston type pressure regulator for a great reliability.

Single stage pressure regulator.

Integrated relief valve protected in the body of the pressure regulator thus ensuring its inviolability and its reliability.

Pressure gauge protected against shocks thanks to a plastic bumper.

1.6 class high pressure gauge (315 bar).

Filter at the inlet to prevent the risks linked to high pressure.

Delivered with a **DEBSON TM2 flowmeter with pre-adjusted flowrates.**

• **Use:**

Pressure regulator :

- Check that the device is undamaged, especially at the inlet connection.
- Slightly and briefly open the gas cylinder knob BEFORE connecting the pressure regulator in order to remove impurities from the cylinder outlet. Then turn it off.

DETENDEURS AVEC DEBITMETRE A ORIFICES PRECALIBRES
PRESSURE REGULATORS WITH FLOWMETER WITH PRE-ADJUSTED FLOWRATES

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
- *Make sure there is a gasket at the inlet of the pressure regulator (for both pin-index and bullnose versions).*
- *Tightly manually screw the bullnose inlet stem on the cylinder knob connection or the pin-index on the cylinder.*
- *Gently open the gas cylinder knob.*

DO NOT USE ANY SPANNER OR PLIERS.

Flowmeter:

- *Screw on a humidifier if necessary, or a nipple to connect the needed tubing for therapy.*
- *Turn the front knob counterclockwise until you reach the requested flowrate.*
- *To stop the flow, turn the gas cylinder knob off.*

Disassembly:

- *First turn the gas cylinder knob off.*
- *Wait until the flow ends or blow off the pressure regulator by opening the flowmeter (the pressure gauge must display 0). Always make sure the pressure regulator is blown off before unscrewing.*
- *Manually unscrew the pressure regulator from the cylinder.*


WHEN UNUSED, DO NOT LET THE DEVICE UNDER PRESSURE.

• **Traceability:**

A unit serial number is engraved on each pressure regulator body for identification and traceability.

8 digits serial number

*1st and 2nd digits = 2 last numbers of manufacturing year
3rd and 4th digits = manufacturing month
3 or 4 last digits = incremented serial number*

Ex : 11 11 1234

• **Cleaning:**

Clean the exterior of the device with water and soap. Rinse and dry. If using disinfecting products, please check their compatibility with plastics. Do not lay under water.

• **Maintenance:**

Device to be serviced every 5 years in our workshop (see next overhaul date affixed to each product).

Only the manufacturer or an authorized person can operate on inner spare parts.

DETENDEURS AVEC DEBITMETRE A ORIFICES PRECALIBRES
PRESSURE REGULATORS WITH FLOWMETER WITH PRE-ADJUSTED FLOWRATES

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• **Safety:**

You are reminded that **the operating personnel must be trained in the handling of gases:**

- Never carry out multiple pressurisations successively.
- Never turn the flowmeter on before opening the High Pressure valve (it should always be set at 0 l/min beforehand).
- Do not turn the flowmeter on directly to maximum flow (gradually increase the flow from 0 to 15 l/min by passing through the intermediary positions).
- Do not close the High Pressure valve with excessive torque (do not force).

Usual precautions with oxygen (O₂):

- **No smoking.**
- **Do not bring into contact with a flame.**
- **Do not lubricate.**
Do not coat the patients' face with greasy substances (Vaseline, ointments).
Handle the equipment with clean hands, free from grease; the wearing of clean examination gloves (latex, nitrile rubber) is recommended.
- Never stand opposite to the valve output when opening but always on the side opposite to the pressure-regulating valve, behind the bottle and taking a step back.
- Never expose the patient to the gas flow.
- Do not use aerosols (hairspray, air-freshener etc.), solvents (alcohol, oil etc.) on the equipment or close by.
- Check for leaks; in the event of a leak, close the valve. Never use a gas cylinder presenting a sealing defect.
- Open the valve gradually.
- Never force the valve when opening.

You are reminded that the safe use of this pressurised medication requires all the information on the label and in the instruction manual to be read carefully and any incident must be declared to the competent authorities.

• **Packing / delivery:**

Delivered in an individual sealed plastic bag with instructions for use.

• **Warranty:**


1 YEAR warranty, parts and service, except in case of damage or accidents due to carelessness, misuse and bad supervision or maintenance.

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	11054	Outlet 9/16"	Brass
	11056	Outlet 1/4G	Brass
	16921	Outlet 1/2"BS	Brass
14	11051	Filter	Bronze
15	11578	Gasket 6.07X1.78	Nitrile
16	11024	Tightening ring	Brass
17	11578	Gasket 6.07X1.78	Nitrile
18	11022	Rotation axis	Brass
19	17715	Tube 01 l/min with stop system	Xenoy
	17717	Tube 05 l/min with stop system	Xenoy
	17714	Tube 15 l/min with stop system	Xenoy
	18525	Tube 25 l/min with stop system	Xenoy
	17716	Tube 50 l/min with stop system	Xenoy
20	17723	Wheel 01 l/min	Xenoy
	17729	Wheel 05 l/min	Xenoy
	17722	Wheel 15 l/min	Xenoy
	18527	Wheel 25 l/min	Xenoy
	17724	Wheel 50 l/min	Xenoy
21	18512	Front disc 01 l/min O ₂	Polypropylene
	18516	Front disc 05 l/min O ₂	Polypropylene
	18511	Front disc 15 l/min O ₂	Polypropylene
	18514	Front disc 25 l/min O ₂	Polypropylene
	18515	Front disc 50 l/min O ₂	Polypropylene
22	11030	Screw	Stainless
23	11031	Gasket	Nylon
24	17811	Spring	Stainless
25	17779	Inox ball	Stainless
26	18478	Disc 01 l/min 4.5 bar	Nickel
	18479	Disc 05 l/min 4.5 bar	Nickel
	18480	Disc 15 l/min 4.5 bar	Nickel
	18482	Disc 25 l/min 4.5 bar	Nickel
	18481	Disc 50 l/min 4.5 bar	Nickel
27	11025	Disc	Brass
28	11030	Screw	Stainless
29	11425	Gasket	Copper
30	11420	Grey bumper for O ₂ pressure gauge	PVC
31	11422	1/8 lateral pressure gauge	Steel