REGULATEUR DE VIDE / VACUUM REGULATOR

## REGULATEUR DE VIDE RVTM3 RVTM3 VACUUM REGULATOR

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Medical device meant to **measure and adjust the vacuum level within the context of medical and surgical suction**s. To be connected to a wall vacuum source. The vacuum regulator should be used with a collection jar and a suction hose.

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• **Classification rule and CE certificate registration No.:**

- Active device of class IIa non sterile - according to the classification of the annex IX of the european directive 93/42/EEC and 2007/47/EEC, chapter 3 rule 11.  
- CE 0197 – TÜV RHEINLAND  
- HD 60020766 0001  
- Directive 93/42/EEC and 2007/47/EEC, Annex II, Article 3

• **UMDNS code and product group:**

13846: "Suction kits"

• **GMDN code and product group:**

36778: "A device that uses negative pressure supplied through the hospital's medical gas supply system, and that is used in the treatment of the patient; e.g. for the aspiration of fluids and/or particulate matter. It is typically used bedside, or in operating theatres."

• **Applied standards to the products:**

**EN ISO 10079: 2009** - Medical suction equipment. Part 3: Suction equipment powered from a vacuum or pressure source.

**EN ISO 7396-1: 2007** - Medical gas pipeline systems - Part 1: Pipelines systems for compressed medical gases and vacuum.

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• **Available suction levels:**

High suction: **From 0 down to - 1000 mbar/hPa**  
(graduated every 50 mbar/hPa)  
**From 0 down to - 760 mmHg/hPa**  
(graduated every 25 mmHg)  
Medium suction: **From 0 down to - 600 mbar/hPa**  
(graduated every 20 mbar/hPa)  
**From 0 down to - 400 mmHg**  
(graduated every 10 mmHg)  
Low suction: **From 0 down to - 250 mbar/hPa**  
(graduated every 10 mbar/hPa)

• **Unit of measurement:**

millibar (mbar) / hectopascal (hPa)  
or millimeter of mercury (mmHg)  
**1 kPa = 10 mbar = 10 hPa = 7.5 mmHg**

• **Reading accuracy:**

Gauge of class 1.6  
Measurement error = 1.6 % of full scale for each value:

250 mbar/hPa = error of +/- 4.00 mbar/hPa  
400 mmHg = error of +/- 6.40 mmHg  
600 mbar/hPa = error of +/- 9.60 mbar/hPa  
1000 mbar/hPa = error of +/- 16.00 mbar/hPa  
760 mmHg = error of +/- 12.16 mmHg

• **Maximum flowrate:**

110 l/min (for a pump flowrate of 160 l/min)

• **Suction limit system:**

0-1000 mbar : without limit  
0- 760 mmHg : without limit  
0- 600 mbar : limit 600 + 40 mbar  
0- 400 mmHg : limit 400 + 30 mmHg  
0- 250 mbar : limit 250 + 40 mbar + mechanical safety valve

• **ON/OFF Switch button:** Push button

• **Regulation knob:**

Central up-front, with a sign for the opening direction  
(counterclockwise)  
Free rotation of the knob at the end of the course (impossible  
blocking).

• **Vacuum gauge:**

Manually adjustable from - 45° to + 45°  
Diameter: 60 mm - Protected by a plastic housing.

• **Connection to the standardized wall outlet:**

Either directly with a direct probe.  
Or with a complete rail mounting system (including rail clamp,  
metal sliding piece, 1.5 m hose and crimped remote probe).

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
• **Outlet connection**

Connect a silicone suction hose of minimum 6.3 mm diameter to the filter up-front.

It is possible to connect the suction hose directly to the safety jar or to the integrated metal tubing nipple of the device in case filters and safety jars are out of stock.

• **Weight with safety jar, without probe:** 390 g.

• **Dimensions with safety jar:**

• **Features:**

Compact and ergonomic device.

Manual adjustment of the vacuum gauge from  $-45^{\circ}$  to  $+45^{\circ}$  for a better visibility.

Quick adjustment of the vacuum level: 2.5 turns are enough to reach the maximum vacuum level.

Free rotation of the regulation knob at the end of the course (impossible blocking).

Supplied with a 100 ml safety jar equipped with a single-use plastic filter up-front. Safety jar made of polycarbonate, autoclavable up to  $134^{\circ}\text{C}$ , unbreakable, equipped with a mechanical anti-overflow system. The safety jar does not require any sterilization except in case of accidental 'liquids' overflow or perforated filter.  
Time and cost savings guaranteed !  
Same filter as for the RVTM2.

Fixing of the safety jar by an easy-click rotation.

Rotation of the safety jar to avoid any pinch of the tubing.

**3 in 1 system – Patented model**

Device equipped with a metal outlet tubing nipple integrated in the body of the regulator. For a better safety emergency suction even be processed if running out of stock of filters or safety jars.

The inner system of the device is protected but accessible at the back for maintenance.

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**Suction air flowrates of vacuum regulators:**

Maximum flowrate of the suction pump: 160 l/min of air

RVTM3 0-250		RVTM3 0-600		RVTM3 0-1000	
Mbar/hPa	l/min	Mbar/hPa	l/min	Mbar/hPa	l/min
25	25	50	75	100	40
50	40	100	100	200	75
100	70	200	110	400	110
150	100	300	110	600	110
200	110	400	110	800	110
250	110	600	110	1000	110

RVTM3 0-400		RVTM3 0-760	
mmHg	l/min	mmHg	l/min
40	75	75	40
75	100	150	75
150	110	300	110
225	110	450	110
300	110	600	110
400	110	750	110

**Packing / delivery:**

Delivered in an individual sealed plastic bag with instructions for use.

**• Adjustment:**

The adjustment of a vacuum regulator must be done in a closed suction network:

- Block the outlet of the vacuum regulator.
- Push the ON/OFF switch button to get the green part visible.
- Turn the regulation knob counterclockwise to increase the suction level to be read on the gauge.
- Get the needle of the gauge to the suction level you wish.

If the adjusted vacuum level is too high, then close the regulation knob and open the vacuum circuit to blow off the device. Then repeat the adjustment operation mentioned above.

**• Closing:**

For a temporary stop, use the ON/OFF switch button to get the red part visible. The vacuum level comes back to its initial value when suction is back on.

For a long-term stop, turn off the regulation knob and the ON/OFF switch button (red part visible). Open the vacuum circuit to blow off the device.

**Visible red part = closed (vacuum is cut)  
Visible green part = open (suction)**

**• Traceability:**

A unit serial number is laser engraved on each vacuum regulator body for identification and traceability.

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*Serial number with 8 digits:*

*1<sup>st</sup> – 2<sup>nd</sup> digits = 2 last numbers of manufacturing year*

*3<sup>rd</sup> – 4<sup>th</sup> digits = manufacturing month*

*Last digits = serial number*

september

Ex: 11 09 0117 → Serial number

2011

• **Cleaning:**

Clean the outside of the device with water and soap. Rinse and dry. Make sure the cleaning products are compatible with plastics. Do not lay under water.

To change the filter and to clean the safety jar, please refer to page 8.

• **Maintenance:**

*Check every 1 to 3 years depending on the use:*

- *Inspection of flow accuracy (vacuum gauge)*
- *Inspection of tightness*

*Once a year, lubricate the ON/OFF switch button with "High Vacuum Grease" silicone grease, ref. 11853. To take the ON-OFF switch button out, pull on the red part.*

*The ON-OFF switch button is protected by a moulding system to avoid any risk of wrong positioning during assembly.*

*To change the vacuum gauge, unscrew the 2 screws at the back of the body.*

*Record every maintenance operation.*

***In case of damage or malfunction, send the device back to the manufacturer or to a registered supplier.***

• **Warranty:**


*1 year warranty, parts and service, except in case of damage or accidents due to carelessness, misuse and bad supervision or maintenance.*

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No	Ref.	Description	Material
<b>VACUUM GAUGE WITH PROTECTION</b>			
1	18675	Back cover for vacuum gauge (white)	ABS
	18581	0-250 mbar vacuum gauge	ABS + brass
	18580	0-600 mbar vacuum gauge	
2	18579	0-1000 mbar vacuum gauge	
	11875	0-400 mmHg vacuum gauge	ABS
	18700	0-760 mmHg vacuum gauge	
3	18674	600/1000 mbar - 400/760 mmHg yellow front cover	
	18743	600/1000 mbar - 400/760 mmHg grey front cover	ABS
	18673	0-250 mbar yellow front cover	
	18742	0-250 mbar grey front cover	
4	15316	Vacuum gauge/nut gasket 7x2.5	Nitrile
5	18676	Vacuum gauge nut	ABS
<b>BODY</b>			
6	11415	Vacuum gauge gasket 13.94x2.62	Nitrile
7	18691	Back cover screw	ABS
8	18669	Yellow back cover	ABS
	18744	Grey back cover	
9	11816	Insert	Brass
10	18694	Valve body	Brass
11	18693	Valve spring	Stainless steel
12	18692	Valve ball	Stainless steel
13	11684	Valve gasket 7.65x1.78	Nitrile
14	18705	Body only	ABS
15	18672	Outlet	Brass
16	18731	Outlet gasket 7.6X2.62	Nitrile
17	18677	Yellow skirt	ABS
	18745	Grey skirt	
18	18686	Yellow front cover	ABS
	18746	Grey front cover	
19	18685	Transparent cover	Polyamide
	18604	Disc 0-1000 mbar	PVC
	18605	Disc 0-600 mbar	
20	18606	Disc 0-250 mbar	
	19490	Disc 0-760 mmHg	
	19491	Disc 0-400 mmHg	
21	11779	Needle gasket 2.9X1.78	Nitrile
22	18683	Needle	Brass
23	18684	Membrane rod	Brass
	18680	Spring 0-250 mbar	Stainless steel
24	18681	Spring 0-600 mbar/0-400 mmHg	
	18682	Spring 0-1000 mbar/0-760 mmHg	
25	11819	Disher	Stainless steel
26	11820	Circlips	Steel
27	18679	Green vacuum regulation knob	ABS
	18756	Grey vacuum regulation knob	
28	18678	Green regulation knob stopper	ABS
	18755	Grey regulation knob stopper	
29	18576	Membrane	Silicone
30	11749	Bumper membrane	Polycarbonate
31	18687	Knob body only	ABS
32	18688	Green pushbutton	ABS
33	18689	Red pushbutton	ABS
34	11839	Knob gasket (x2) 6.75x1.78	Nitrile
42	18789	Gasket for switch-button	Polyurethane
	11826	Inlet connector 12x100 F	Brass
43	11823	Inlet connector 1/4 G M	Brass
	11825	Inlet connector 1/8 NPT F	Brass
<b>100 ml SAFETY JAR</b>			
35	11780	Safety jar gasket 44.04X3.53	Nitrile
36	18690	Cover only	Polypropylene
37	17294	Safety cage	Polypropylene
38	11698	Overflow ball	Polypropylene
39	17630	Safety jar only	Polycarbonate
40	11701	Safety gasket 9.19X2.62	Nitrile
41	11813	Tube of 10 single-use antibacterial plastic filters	Borosilicate fibreglass + laminated polyester + polypropylene

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• **Capacity: 100 ml**

• **Mechanical anti-overflow safety system**

• **Fixing and removing of the safety jar:**

Fixing of the safety jar under the vacuum regulator by an easy-click rotation:

Line up the triangular mark of the safety jar lid to the one on the skirt of the vacuum regulator. Then turn the safety jar of one quarter to the right.

Removing:

Turn the safety jar to the left. Line up the two triangular marks and pull down the safety jar.

• **Cleaning of the safety jar:**

The safety jar being protected against infections by the filter up-front, it is not necessary to disinfect it except in case of accidental liquids' overflow or perforated filter (see below for disinfection procedure).

Disinfection procedure:

- Unscrew the safety jar and disassemble the cover from the jar.
- Dump the bottle in a decontaminating and cleaning solution.
- Clean, rinse thoroughly, let soak for 1 minute in clean water and dry.
- Sterilize or disinfect.
- Put a new filter at the right place.
- Re-assemble the parts together and re-screw the safety jar on the vacuum regulator.

**The safety jar withstands the autoclave up to 134°C.**

Caution: the compatibility of the cleaning products used with polycarbonate must be checked. For information, products containing amphoteric components, glutaraldehyde, formaldehyde, glyoxal, succinic acids, **are not compatible** with polycarbonates. Also avoid solutions for floors and surfaces.

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• **Single-use antibacterial plastic filter up-front**

Non sterile filter  
Product class I

Dimensions:  
Diameter: 50 mm  
Thickness: 0.47 mm

Performance:  
Density: 100 g/m<sup>2</sup>  
DEHS retention: 5.3 cm/s  
Particle diameter: 0.3 µm  
Efficiency: > 99.97 %

Material:  
Filter: borosilicate fiberglass + laminated polyester  
Plastic housing: polypropylene

Packing: Tube of 10 single-use antibacterial plastic filters

**Change of the filter:**

Change the antibacterial filter after each patient or, for a long-stay patient, after inspection of the filter clogging level.

**Fixing of the filter:**  
Push until click-lock adjustment

**Removing of the filter:**  
Tightly pull the filter while operating a rotational movement.

Throw the filter away respecting the usual precautions.