VENTURI TM / VENTURI TM EJECTOR

VENTURI TM VENTURI TM EJECTOR

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Medical device meant to adjust and measure a suction from a source of pressured gas on the wall or from a cylinder. It enables to suck in liquids or mucus in the absence of vacuum pipeline network. It should be connected to a source of Oxygen or Medical Air. The VENTURI ejector should be associated with a collection jar and a suction hose.

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• **Classification rule and CE certificate registration No.:**

*- Active device of class IIa non sterile - according to the classification of the annexe IX of the european directive 93/42/EEC and 2007/47/EEC, chapter 3 rule 11.
- CE 0197 – TÜV RHEINLAND
- HD 60004355 0001
- Directive 93/42/EEC and 2007/47/EEC, Annexe II, Article 3*

• **UMDNS code and product group:**

13846: "Suction kits"

• **GMDN code and product group:**

36778: "A device that uses negative pressure supplied through the hospital's medical gas supply system, and that is used in the treatment of the patient; e.g. for the aspiration of fluids and/or particulate matter. It is typically used bedside, or in operating theatres."

• **Applied standards to the products:**

EN ISO 10079: 2009 - Medical suction equipment.
Part 3: Suction equipment powered from a vacuum or pressure source.

EN ISO 7396-1: 2007 - Medical gas pipeline systems - Part 1: Pipelines systems for compressed medical gases and vacuum.

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
• **Available gases**: Oxygen (O₂) and Medical Air

• **Working pressure**: 4.5 bar ± 0.5

• **Unit of measurement**: millibar (mbar)

• **Reading accuracy**:

Gauge of class 1.6 (1.6 % of measurement error).

• **Suction flow**: 25 l/min at - 650 mbar

• **Consumption**: 10 % in l/min of the adjusted vacuum value
(Ex: 200 mbar → consumption of 20 l/min)

• **Maximum suction level**: - 650 mbar ± 10 %

• **Working temperature**: 23°C (73.4°F)

• **Connection to the wall outlet**:

Using a direct probe or a rail mounting system.

• **Weight (without probe)**: 650 g

• **Dimensions (without probe)**:

• **Technical data**:

Square and wide vacuum gauge, graduated in mbar
(0-1000 mbar), for an easy reading.

Gauge: 48 mm x 48 mm.

Regulation knob: Needle valve knob.

Anodized aluminium block ensuring a great robustness.

Ejector calibration cannot be disturbed.

Supplied in standard with a 100 ml safety jar with plastic filter up-front to protect the device against any liquids' overflow. Also available with:

- A 150 ml or 500 ml safety jar with paper filter in the cover.
- An outlet tubing nipple.

Safety jar made of polycarbonate, autoclavable up to 134°C, unbreakable and fitted with an anti-overflow safety device.

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• **Use:**

*The adjustment of the VENTURI ejector must be carried out in closed circuit.
Block the outlet of the VENTURI.
Gradually turn the regulation knob to the left till the needle indicates the requested suction level on the vacuum gauge.*

• **Traceability:**

A unit serial number is engraved on each product body for identification and traceability.

8 digits serial number:

- 1st and 2nd digits = 2 last numbers of manufacturing year
- 3rd and 4th digits = manufacturing month
- 3 or 4 last digits = incremented serial number

• **Cleaning:**

*Clean the outside of the device with water and soap. Wipe and dry. If using disinfecting products, check their compatibility with plastics (polypropylene and polycarbonate).
Do not use cleaning products for surfaces.*

To change the filter and to clean the safety jar, please refer to page 5.

• **Maintenance:**

Check every 1 to 3 years depending on the use:

- inspection of flow accuracy (gauge),
- inspection of tightness.

Replace the silent filter once a year.

Change the antibacterial filter of the safety jar after each patient.

• **Safety:**

*Make sure to ventilate carefully the room to avoid hyper-oxygenation.
Do not use close to any ignition point.
Do not bring the device into contact with greasy products.*

• **Packing / delivery:**

Delivered in an individual sealed plastic bag with instructions for use.

• **Warranty:**


1 year warranty, parts and service, except in case of damage or accidents due to carelessness, misuse and bad supervision or maintenance.

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- **Capacity**: 100 ml as standard
150 ml on request

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- **Mechanical overflow safety system**

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Current model: Single-use antibacterial plastic filter up-front.

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Former model: Single-use antibacterial paper filter in the cover.

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See on page 8 the technical data of the antibacterial filter.

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1/ Antibacterial plastic filter up-front:

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Change the filter after each patient or, in case of a long stay patient, when necessary after inspection of the filter clogging level. Pull the filter with a rotating movement and throw it away respecting the safety rules. Put a new filter by pushing it until ratchet mechanism.

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The safety jar being protected against infections by the filter up-front, it is not necessary to disinfect it any longer except in case of occasional liquids' overflow or perforated filter (see below for disinfection procedure).

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2/ Antibacterial paper filter in the cover:

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Change the filter after each patient or, in case of long stay patient, when necessary, after inspection of the filter clogging level. Each time you change the filter it is essential to disinfect the safety jar according to the procedure below:

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- Unscrew the safety jar and disassemble the cover from the jar.
- Dump the bottle in a decontaminating and cleaning solution.
- Clean, rinse thoroughly, let soak for 1 minute in clean water and dry.
- Sterilize or disinfect.
- Put a new filter at the right place.
- Re-assemble the parts together and re-screw the safety jar on the vacuum regulator.

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The safety jar withstands the autoclave up to 134°C.

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
Caution: the compatibility of the cleaning products used with polycarbonate must be checked. For information, products containing amphoteric components, glutaraldehyde, formaldehyde, glyoxal, succinic acids, **are not compatible** with polycarbonates. Also avoid products for floors and surfaces.

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• **Technical view**: Safety jar with antibacterial plastic filter up-front

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
No	Ref.	Description	Materials
1	11772	COMPLETE COVER WITH ANTIBACTERIAL PLASTIC FILTER UP-FRONT	
1A	11781	1/2G inlet gasket	Ethylene Propylene Diene Monomer
1B	11845	1/2G nut	Brass
1C	11844	Inlet stem	Brass
1D	11770	Cover only with filter holder	Polypropylene
1E	11701	O-ring for safety cage	Ethylene Propylene Diene Monomer
1F	11780	O-ring for safety jar	Ethylene Propylene Diene Monomer
1G	11698	Overflow safety ball	Polypropylene
1H	11697	Overflow safety cage	Polypropylene
JARS ONLY			
2	17630	100 ml jar only	Polycarbonate
	11731	150 ml jar only	Polycarbonate

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• **Technical view**: Safety jar with antibacterial paper filter in the cover

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No	Ref.	Description	Material
1	11741	COMPLETE COVER WITH ANTIBACTERIAL PAPER FILTER IN THE COVER	
1A	11781	1/2G inlet gasket	Ethylene Propylene Diene Monomer
1B	11845	1/2G nut	Brass
1C	11844	Inlet stem	Brass
1D	11743	Filter translucent cover only	Polycarbonate
1E	11780	O-ring for safety jar	Ethylene Propylene Diene Monomer
1F	11817	Antibacterial paper filter	Borosilicate fibreglass + laminated polyester
1G	11744	Safety cover only	Polycarbonate
1H	11811	Plunger tube	Polycarbonate
1I	11701	O'ring for safety cage	Polypropylene Ethylene Propylene Diene Monomer
1J	11698	Overflow safety ball	Polypropylene
1K	11697	Overflow safety cage	Polypropylene
JARS ONLY			
2	17630	100 ml jar only	Polycarbonate
	11731	150 ml jar only	Polycarbonate

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SINGLE-USE ANTIBACTERIAL FILTERS FOR SAFETY JARS

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• **Material:**

- Filter: borosilicate fibreglass + laminated polyester
- Plastic housing: polypropylene

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Ref.: 11813
Tube of 10 plastic filters
for current safety jar with
filter up-front

Ref.: 11818
Box of 100 paper filters for
former safety jar