# RV™ 3

# CE 0459

FR	Notice d'instructions
EN	Instructions for use
DE	Bedienungsanleitung
ES	Folleto de instrucciones
BG	Инструкции за експлоатация
EL	Οδηγίες χρήσης
HU	Használati útmutató
LT	Naudojimo instrukcija
MT	Struzzjonijiet għall-użu
PL	Instrukcja obsługi
RO	Manual de Instrucțiuni
SK	Návod na použitie
SL	Navodila za uporabo
CS	Návod k použití
DA	Brugsanvisning
AR	إرشادات الاستخدام
IT	Avvertenze e istruzioni
NL	Gebruiksaanwijzing
PT	Manual de Instruções
SV	Bruksanvisning

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1. Presentation
2. Versions
3. Technical data
4. Instructions for use
5. Symbols
6. Safety advice
7. Cleaning and disinfection
8. Transport and storage
9. Maintenance
10. Guarantee

## 1. Presentation

This medical device **measures and adjusts the vacuum level during medical or surgical aspiration procedures**. It connects to a wall-mounted vacuum outlet. The vacuum regulator must be used in conjunction with a collection jar and a suction tube.

User profile:

Users are medical staff: doctors, nurses.

## 2. Versions

The RVTM3 vacuum regulator is available in the following versions:

- Maximum vacuum levels available: 0-250, 0-600, 0-1000 mbar/hPa and 0-400, 0-760 mmHg.
- Connector: compliant with NF S 90-116 - BS 5682 - DIN 13260 - US Ohmeda Diamond - NORDIC - UNI.
- Direct connection or via a rail mounting (tube + crimped tube holder + rail clamp).
- Outlet: tubing nipple or tubing nipple + safety jar
- Available in single or double versions.

## 3. Technical data

- Complies with standard ISO 10079-3
- Units: millibar (mbar), hectoPascal (hPa) and millimetre of mercury (mmHg): 10 mbar = 10hPa = 7.5 mmHg.
- Measurement accuracy: 1.6 % of the full scale indicated on the vacuum gauge.
- Vacuum gauge can swivel through 90° (from - 45° to + 45°).
- The regulation knob rotates freely when it reaches the end of its control range (it cannot be excessively tightened against a stop).
- Safety jar can swivel through 90° (from - 45° to + 45°).
- Antibacterial filter on front of device: 99.97% efficiency for 0.3 µm particles.
- Individual serial number: laser-engraved on the device's body.
- Safety valve fitted for the 0-250 mbar version.
- Metal connection nipple built into the body of the vacuum regulator.
- Maximum flow rate: 110 l/min (when used with a pump generating a flow rate of 160 l/min).
- Calibrated admission valve (3 l/min when in service) to allow dynamic adjustment of the vacuum level.
- Ambient temperature for use: 23°C.

## 4. Instructions for use

- The suction equipment must only be used by suitably-trained personnel.
- Influence of the source of vacuum: caution, check the maximum vacuum level which the source of vacuum can produce. The maximum vacuum level set by the vacuum regulator cannot be greater than the maximum value produced by the source.

### **Vacuum regulator**

- Make sure that the device is undamaged and that the connector is compatible with the source connection.
- Make sure that the regulation knob is in its closed or "zero-suction" position.
- Switch the device to its OFF position (indicated by the red button protruding from the device).
- Connect the device to the source of vacuum outlet.
- Connect the vacuum regulator's outlet (tubing nipple or safety jar) to the collection jar (Ø 9 mm) using a suction tube with a minimum diameter of 6.3 mm.
- To prepare for aspirating, press on the red button (so that the green button protrudes).
- By turning the regulation knob on the front of the device anti-clockwise, set the vacuum level with the patient circuit closed off.

➔ To interrupt the aspiration temporarily, use the ON/OFF button. The vacuum level returns to the value set previously when the source is restored.

➔ To switch off the device for a long period, set the regulation knob to zero flow and ensure that the OFF (red) button is protruding (so that the circuit can drain).

### **Safety jar**

- **Fitting:** to fit the jar under the vacuum regulator: line up the mark on the cover of the safety jar with the mark on the skirt of the body of the vacuum regulator then rotate the safety jar to the right through ninety degrees.
- **Removing:** rotate the safety jar to the left until the mark on the cover of the safety jar lines up with the mark on the skirt of the body of the vacuum regulator then pull the jar downwards.


#### **Why use an antibacterial filter?**

An aerosol may be produced during suction which might contaminate the vacuum regulator, connectors, wall outlets, lines and vacuum pumps. Moreover, if there is no flow, and no filter present, bacteria could move freely through the patient circuit.

### **Filter**

- **Fitting:** push in fully until it clicks into place
- **Removing:** pull and twist at the same time.
- Dispose of the filter by observing the appropriate precautions.

## 5. Symbols

	Manufacturer		Conformity with European directive 93/42/EEC, notified body No.0459	
	Refer to the instructions		Refer to the accompanying documents	 Adjustment of the vacuum level
	Reference No. on packaging		Lot No. on packaging	 Single use
	Serial No. of the product, with the format YYMMxxxx (YY: year of manufacture, MM: month of manufacture, xxxx: unique number)			Use-by date, with the format YYYY/MM (YYYY: year, MM: month)

EN

## 6. Safety advice

- Always test the operation of the device before using it on the patient.
- The vacuum gauge with a dial display is sensitive to impacts. After an impact or if it is dropped, check that the value displayed is correct.
- Make sure that the needle indicates 0 when the device is not being used.
- Make sure that a filter is always fitted.
- Always disconnect the vacuum regulator at the connector if liquid or solid matter has been aspirated into the vacuum pump.
- Caution: the length of the tubing may have an impact on suction performance.
- It is not advisable to use the vacuum regulator in an MRI environment.

## 7. Transport and storage

- To prevent zero-offset errors affecting the vacuum gauge's needle, the device should be laid horizontal for transport.
- Store at between -20 and +60°C in a clean and dry place.
- Keep the instructions with the packaging for as long as possible.

## 8. Cleaning and disinfection

Under normal conditions of use, the RVTM3 does not require disinfection since the safety jar is protected by the filter on the front of the device. **Change the filter after each patient.**

### Vacuum regulator

Use soapy water to clean the outside of the device. Rinse thoroughly and dry. Do not immerse in water or other liquid.

### **Caution:**

- When using decontaminant products, make sure they are compatible with plastic.
- Do not use surface cleaners.

### Safety jar

**Since the safety jar is protected by an upstream filter on the front of the device, there is no longer any need to sterilize it: simply replace the filter after each patient or, if the patient is admitted for a long period, replace as necessary after inspecting the filter to check the degree of fouling.**

**However, if the liquid accidentally overflows or if the filter perforates, then the safety jar must be disinfected or sterilized as follows:**

1. Remove the safety jar and disassemble the cover from the jar.
2. Immerse safety jar parts in a pre-disinfecting solution (respecting the manufacturer's instructions concerning soaking, rinsing and drying times).
3. Clean, rinse thoroughly, soak for 1 minute in clean water and then dry.
4. Sterilize or disinfect in compliance with the instructions provided by the manufacturer of the particular product used.
5. Reassemble each part and refit the safety jar.

**The RVTM3 safety jar can withstand at least 30 cycles of autoclaving at 134°C.**

## 9. Maintenance

No preventive maintenance is required on the RVTM3 during its service life.

Technologie Médicale recommends inspecting the device every 1 to 3 years, depending on how often it is used:

- Check the vacuum level displayed (on the vacuum gauge).
- Check its operation.

### Spare parts:

- Safety jar
- Tube of 10 filters


Contact the distributor to obtain the product reference codes corresponding to the vacuum regulator used.

**Only use TECHNOLOGIE MEDICALE spare parts.**

**In the event of damage or malfunction, return the device to the manufacturer or to an approved repairer.**

## 10. Guarantee

- The device is covered by a 1-year warranty (parts and labour), which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.
- Service life: 10 years.
- Date when **CE** marking first affixed: 2011

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