BLENDER TM2

CE 0459

FR	Notice d'instructions	p. 2
EN	Instructions for use	p. 4
DE	Bedienungsanleitung	p. 6
ES	Manual de instrucciones	p. 8
PT	Manual de Instruções	p. 10
IT	Avvertenze e istruzioni	p. 12
NL	Gebruiksaanwijzing	p. 14
SV	Bruksanvisning	p. 16
AR	إرشادات الاستخدام	p. 18
BG	Указание за употреба	p. 20
EL	Εγχειρίδιο οδηγιών	p. 22
HU	Útmutató	p. 24
LT	Naudojimo instrukcijos	p. 26
MT	Manwal tal-istruzzjonijiet	p. 28
PL	Instrukcja użytkowania	p. 30
RO	Instrucțiuni de utilizare	p. 32
SK	Návod na použitie	p. 34
SL	Navodila za uporabo	p. 36
CZ	Návod k použití	p. 38
DA	Brugsanvisning	p. 40
FI	Käyttöohjeet	p. 42
LV	Lietošanas instrukcija	p. 44
NO	Bruksanvisning	p. 46
TR	Kullanım talimatları	p. 48

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1. Presentation

The BLENDER TM2 is a device that mixes air and oxygen. It has been designed to dispense a mixture of medical air and oxygen to newborn babies, children and adults. It is intended to be used by qualified and trained personnel, under a doctor's supervision, in a professional healthcare setting, e.g. a hospital. The BLENDER TM2 lets the doctor select the oxygen concentration to be dispensed to the patient and the flow rate of this mixture.

The BLENDER TM2 has a floating ball flowmeter (rotameter) permanently fitted to the left-hand and/or right-hand side of the device to ensure that the perfect concentration of gases is delivered at low flow rates (RTM3 Fuite). The RTM3 Fuite has a 3 or 10 L/min bleed feature that is activated when the flowmeter is opened for use and disabled when the flowmeter is closed off; in all other respects this flowmeter operates in the same way as a normal RTM3.

The purpose of these instructions for use is to provide doctors, healthcare personnel and patients/users with the information needed to ensure that the device is used correctly.

Please read the instructions for use carefully before using the device for the first time!

2. Versions

The following versions of the BLENDER TM2 are available:

- Versions: Low Flow / High-Low Flow / High Flow
 - Gas: AIR-OXYGEN mixture
 - Pressure at device inlet: O₂ 4.5 bar (+/- 0.5 bar) Air 4.5 bar (+/- 0.5 bar)
 - Inlet connector: NIST
 - Connection via a rail mounting
 - Outlet connector: Male DISS / Flowmeter outlet / Olive / tubing nipple
 - Flow ranges available:
 - LOW FLOW version with RTM3 Fuite: 0-5 l/min; 0-15 l/min; 0-30 l/min LOW FLOW version with RTM3: 3-30 l/min
 - LOW FLOW version with RTM3: 5-30 l/min
 - HIGH-LOW FLOW version with RTM3 Fuite: 2-15 l/min; 2-30 l/min; 2-70 l/min HIGH-LOW FLOW version with RTM3: 15-70 l/min
 - HIGH FLOW version with RTM3: 15-70 l/min
 - FiO₂ accuracy for the LOW FLOW versions:
 - +/-3% with bleed flow for a flow rate of 3 l/min or higher
 - +/-3% without bleed flow for a flow rate of 6 l/min or higher
- Outside these flow ranges and for other versions, the Blender TM2 is accurate to +/-5% in compliance with ISO 11195 requirements


3. Technical data

The BLENDER TM2 has an aluminium body, protective cover and control knob and a polyethylene front face.

- Complies with standard EN ISO 11195
- Oxygen % range: 21 to 100 %
- Pressure at device inlet: 4.5 bar +/- 0.5
- The performance of the BLENDER TM2 is unreliable if either inlet pressure is less than 4 bar or more than 5 bar. Do not connect the device to a patient if either inlet pressure is outside this range.
- The accuracy of the indicated values conforms to NF EN ISO 11195.
- Maximum flow rate for Low Flow version: ≥30 l/min (refer to the table of technical data).
- Maximum flow rate for High-Low Flow version: ≥120 l/min (refer to table of technical data).
- Alarm/Bypass Reset: when the pressure difference between the inlet gases is ≥ 2 bar.
- Inlet filters: mesh size of 48 µm.
- Serial number: with year, month of manufacture and individual number.
- Pressure-compensated flowmeter
- Accuracy of the indicated values for the RTM3 Fuite: ± 10 % of the indicated value or ± 0.5 l/min, whichever is the highest.
- Ambient temperature for use: 23°C.
- Backflow of gas ≤ 10 ml/h (Clause 9 of standard ISO 11195; Test described in 4.2)
- A third outlet provided under the device can be used to connect an O₂ analyser.

Model	Flow rates with no bleed	Flow rates with bleed
Low Flow	5 to 30 l/min	0-30 l/min (bleed of 3 l/min)
High-Low Flow	15-120 l/min (no bleed)	2-108 l/min (bleed of 10-12 l/min)
High Flow	15-120 min	/

The flow with no bleed outlet from the Blender TM2 is indicated by a label on the device (refer to 5. Symbols).

4. Instructions

The following checks must be made by a qualified technician before the BLENDER TM2 is used for the first time.

NB: Do not use the blender if it fails these tests; call the manufacturer's after-sales service.

Test the performance of the device every time before use and/or if the device is dropped (refer to 4.3).

4.1 Start by connecting the 4.5-bar air and oxygen supplies to the device (Air INLET and O₂ INLET) and set the blender to 60% (the alarm should not be triggered). Check that the oxygen concentration is indeed 60% using a calibrated oxygen analyser that conforms to ISO 80601-2-55. Disconnect the oxygen supply to the Blender TM2 and check that the alarm is triggered. Once it has been triggered, reconnect the oxygen supply and check that the alarm stops. Re-check the oxygen concentration. Disconnect the air supply to the Blender TM2 and check that the alarm is triggered again. Once it has been triggered, reconnect the air supply and re-check the oxygen concentration. The device's alarm is also triggered in the event of an abnormal change in gas pressure of more than 1.5 bar at either of the inlets.

4.2 To check the backflow of gas, connect a tube to the device's air inlet. Place the end of the tube in a glass of water and supply the device with oxygen. No bubbles should be seen in the glass. Repeat this operation on the oxygen inlet.

4.3 Connect the gas outlet from the Blender TM2 (Air/O₂ OUTLET), directly or via a high-pressure tube to the equipment which the blender is to be used with. Set the control knob to the required oxygen concentration. Check the inlet pressures (4.5 bar) and adjust the equipment used. Use a calibrated oxygen analyser to check that the gas delivered to the patient corresponds to the value indicated by the control knob. If the oxygen concentration has changed, wait for 60 seconds (to give the system time to balance) before repeating the check using the analyser.

If using a Blender TM2 High-Low Flow at a low flow rate, use the outlet fitted with the RTM3 Fuite. Set the control knob to the required FiO₂. Open the Air and O₂ inlets, and then the RTM3 Fuite. Check that the indicated concentration is correct using a calibrated oxygen analyser.

5. Warnings

All BLENDER TM2 devices are tested in our manufacturing plant. The tests performed are presented in the appended validation table for the BLENDER TM2.

The BLENDER TM2's alarm is triggered by any abnormal change in inlet gas pressure of more than 1.5 bar. The activation of the alarm will modify both the outlet flow rate from the blender and the concentration of gases.

The BLENDER TM2's alarm is also triggered when the supply of one of the device's two inlet gases is interrupted.

The BLENDER TM2 is designed to operate in the supply pressure ranges defined in the technical data. The conditions specific to the outlet with no bleed are indicated on the device by the label shown in the Symbols section.

The BLENDER TM2 is not compatible with the use of Oxygen 93

The BLENDER TM2 was degreased prior to delivery so that it can be used with oxygen.

Never ventilate a patient if no monitoring or remote monitoring is in place.

The hoses supplied (if requested by the customer) with the BLENDER TM2 conform to standard ISO 5359.

Moisture or contamination may affect the operation of the blender. The gas supplied must be always be clean and dry. The air supplied must conform to the USP standard for compressed air and to standard ANSI Z86.1-1973 F and the dew point of the water vapour must not be more than 5° F (- 15° C) to prevent condensation forming at the lowest ambient temperature to which the blender and its accessories are exposed. The oxygen supply must meet the requirements for medical-grade oxygen stated in the USP. The two gases must contain < 37.5 milligrams of water per cubic meter of gas.

6. Symbols

	Do not use grease, lubricant or oil		Refer to the instructions
	Manufacturer	YYMMxxxx	Serial No. engraved on the device: YY: Year of manufacture, MM: Month of manufacture, xxxx: Individual number
Revision mm-yyyy checking	Mention of the preventive maintenance service date recommended by Technologie Médicale mm: month in which service is due, yyyy: year in which service is due	MIN ≥ x L/MIN MAX ≥ x L/MIN	Indicates the flow with no bleed outlet on the device MIN = Minimum MAX = Maximum X = flow rate value
	Made and assembled in France		Indicates the RTM3 Fuite with an outlet of x L/min X = bleed flow rate.
x - x L/MIN	Indicates the outlet fitted with an RTM3 Fuite. To be used in the flow rate range from x to x L/min X = flow rate		Disposal of the equipment should comply with hospital recommendations

7. Safety advice

- To prevent any risk of fire or explosion:
 - ✓ Do not smoke near the device.
 - ✓ Do not operate the device near a source of ignition.
 - ✓ Do not let the device come into contact with greasy or oily products.
 - ✓ Do not handle the device with greasy or oily hands.
 - ✓ Do not disassemble if the device is pressurised.
- Do not create an obstruction to restrict the outlet from the medical device.
- Effect of pressure: Caution: check the inlet pressure fed to the device. If this inlet pressure is not within the range indicated on the device, there is a danger that the flow rate delivered might be outside the permitted normative range.
- Effect of temperature: Caution: the accuracy of the flow rates delivered by the device varies as a function of the ambient temperature (between 0°C and 40°C).
- Do not disassemble the device's inlet or outlet connectors.
- Ensure that the device is always in a vertical position.
- Caution! Do not allow water from a connected humidifier (or other item of equipment) to enter the device when it is transported; always unscrew the humidifier before handling the device in any way.
- Do not obstruct the audible alarm outlet indicated under the device.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- When using lubricant products, make sure they are compatible with oxygen.

8. Cleaning and disinfection

- Disconnect the device from the gas supplies. Clean the outside of the device with a clean cloth, dampened with water. Leave the device to dry. When using decontaminant cleaning products, make sure they are compatible with the plastic used to make the front face of the device (polyethylene). Make sure that nothing is inserted into the device's alarm orifice.
- Do not immerse in any liquid.
- Do not use surface cleaners

9. Storage

Store at between -20 and +40°C in a clean and dry place. Store in its original packaging, if possible.

Humidity and contamination may affect the operation of the blender.

10. Maintenance

Only individuals authorised by Technologie Médicale may disassemble the device. A preventive maintenance service should be performed every two years.

Perform an annual inspection of the device (in-between every preventive maintenance service):

- ✓ Check the device's inlet connections.
- ✓ Check that the gas is delivered correctly and at the right concentration at a range of flow rates.
- ✓ Check that there are no leaks when the flowmeter's valve is closed.
- ✓ Perform the tests indicated in 4. Instructions for use.

11. Guarantee

The device is covered by a 1-year warranty (parts and labour), which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance. Every 2 years, a preventive maintenance service should be carried out to ensure that the device can continue to operate permanently and satisfactorily. This frequency is determined by several factors, such as the frequency and duration of use, the quality of the gases handled and the environmental conditions.

All repair work must be performed by Technologie Médicale or by an authorised technical agent. Technologie Médicale disclaims all liability for the consequences of any unauthorised technical work on any BLENDER TM2.

The date recommended by Technologie Médicale for the preventive maintenance service is shown on the label (refer to 5. Symbols).

The preventive maintenance service, carried out in Technologie Médicale's workshops, includes the following work: General servicing, replacement of parts susceptible to wear, cleaning in an ultrasonic tank and calibration of the product.

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Date d'apposition du premier marquage © 2021