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Collection jar, reusable

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1. Presentation

This medical device is used to collect secretions and other substances when performing medical section on a patient. The device must be connected, via a suction tube, to a fully-equipped suction system incorporating an electrical or mechanical suction unit.

User profile:

Users are medical staff: doctors, nurses.

2. Versions

The collection jar is available with the following options:

- Capacities: 1 L, 2 L or 4 L.
- Materials: polysulfone or polycarbonate.
- Lid: screw-on metal lid, press-on plastic lid, press-on silicone ring with metal central plate, screw-on plastic lid.

3. Technical data

- Complies with EN ISO 10079-3
- Units: millilitre (ml).
- Non-return safety system built into the lid.
- Integrated hanger for the 2 L and 1 L short jars.
- Reusable.
- Ambient temperatures for use: 0 to 40°C.

4. Instructions for use

- Check that the lid is firmly and correctly seated on the jar.
- Check that the "patient" and "vacuum" tubing nipples are tightened correctly.
- Connect up the jar using vacuum tubing (\varnothing 6 to 8 mm for 1 and 2 L jars. - \varnothing 8 to 10 mm for 4 L jars):
 - o On the patient side: to the suction probe.
 - o On the vacuum side: to the vacuum regulator or electrical suction unit.
- Check that each component of the assembly is sealed (before fitting the press-on plastic lid, check that the base of the lid has been inserted correctly into the groove in the seal; this groove may need to be eased open so that the lid can be fitted correctly into the groove)

When the suction is no longer required, do not handle the jar and or its lid until you have checked that the vacuum source is disconnected and that there is no pressurisation of the circuit.

The collection jar is fitted with an overflow protection device, consisting of a ball designed to float up into an orifice to seal it. When this overflow device is activated the suction is stopped. If this should occur, disconnect the vacuum source and then disconnect all tubing or probes from the jar. Remove the lid and empty the jar into a suitable disposal location. Clean and disinfect the jar as described opposite.

5. Symbols

	Refer to the instructions
	Manufacturer

6. Safety advice

- This device must be handled and used by trained and qualified personnel.
- Maintain the jar in a vertical position when in use.
- Always test the operation of the device before using it on the patient.
- Check that a filter is fitted at some point in the vacuum circuit (collection jar or vacuum regulator).
- Do not handle the device when it is depressurised.
- The collection jar is not compatible with use in an MRI environment.
- If liquid flows past the jar into the suction system, withdraw the jar from service and arrange for it to be checked.

7. Cleaning and disinfection

1. Disassemble the components of the jar and lid
2. Immerse the parts in a decontamination or pre-disinfecting solution (respecting the manufacturer's instructions concerning dilutions and soaking, rinsing and drying times).
3. Clean (do not use abrasive products or scourers), rinse thoroughly, soak for 1 minute in clean water and then dry
4. Sterilise or disinfect in compliance with the manufacturer's instructions (keep parts separate, keep jar vertical with its base at the top, do not stack equipment on top of the jars)
5. Reassemble the device

Caution, the rigidity of the jars is impaired by repeated washing in automatic cleaning equipment.

Caution, check that the products used are compatible with polycarbonate or polysulfone.

**Polysulfone can be autoclaved at temperatures of up to 134°C.
Polycarbonate can be autoclaved at temperatures of up to 121°C.**

8. Storage

- Store at between -20 and +60°C in a clean and dry place
- Keep the packaging for as long as possible

9. Maintenance

Every time before the device is used, check that the jar assembly seals correctly and check the condition of the seals.

Spare parts	<u>Use exclusively</u>
Contact us for details about these parts.	TECHNOLOGIE MEDICALE spare parts

10. Guarantee

The device is covered by a 1-year warranty (parts and labour), which does not cover normal wear and tear, damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.

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Datum prvního označení © 2000