CE 0459

DEBFLO DEBPLUS

FR	Notice d'instructions
EN	Instructions for use
DE	Bedienungsanleitung
ES	Folleto de instrucciones
PT	Manual de Instruções
IT	Avvertenze e istruzioni
SV	Bruksanvisning
NL	Gebruiksaanwijzing
AR	إرشادات الاستخدام
BG	Инструкции за експлоатация
CS	Návod k použití
DA	Bruksanvisning
EL	Οδηγίες χρήσης
FI	Käyttöohjeet
HU	Használati útmutató
LT	Naudojimo instrukcija
MT	Struzzjonijiet għall-użu
PL	Instrukcja obsługi
RO	Manual de Instrucțiuni
SK	Návod na použitie

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1. Introduction

DEBFLO and DEBPLUS are medical devices that regulate the flow of a gas. They connect to a wall-mounted source of pressurised gas or to a rapid outlet of a pressure regulator.

User profile:

Users are medical staff: doctors, nurses.

2. Versions

DEBFLO and DEBPLUS are available in the following versions:

- Gas: Oxygen (O₂) or medical air.
- Inlet connection: in compliance with the following standard: NF S 90-116 - BS 5682 - DIN 13260 US Ohmeda Diamond - NORDIC - UNI, etc.
- Direct connection or rail mounted (hose + crimped hose holder end cap + rail clamp).
- Threaded outlet: 12x125 M - 9/16" M - 1/4G M - 1/2"BS F.
- Single or double version.
- Flow ranges in l/min (15 flow positions):

0 - 5 l/min	0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	1	1.5	2	2.5	3	4	5
0 - 15 l/min	0	0.5	1	1.5	2	2.5	3	3.5	4	5	6	7	8	9	12	15
0 - 50 l/min	0	3	5	6	8	9	10	12	15	20	25	30	35	40	45	50

3. Technical characteristics

- Compliant with the standard NF EN ISO 15002 for DEBFLO.
- Compliant with the standard NF EN ISO 10524-4 for DEBPLUS.
- Class IIa medical device.
- Unit of measurement: litres per minute (l/min).
- Supply pressure: 4.5 bar ± 0.5 bar (for DEBFLO); 2.8 to 8 bar (for DEBPLUS).
- Integrated pressure regulation for DEBPLUS.
- DEBPLUS outlet pressure: 2.2 bar max.
- Accuracy of the measured values: ± 30% ≤ 1.5 l/min; ± 20% > 1.5 l/min.
- Continuous flow between two positions.
- Inlet filter: pore size of 35 µm.
- Unit serial number with year and month of manufacture.
- Operating temperature: 0°C to 40°C.

4. Method of use

- Check whether the device is in good condition.
- Ensure that the flow rate is set to 0.
- Connect the device to the wall-mounted outlet or to the rapid outlet of a pressure regulator.
- Screw in a humidifier if required, or a nipple so that the oxygen tubing can be connected to it.
- Turn the regulation knob on the front side to the left until the desired flow rate is displayed.
- Attach the tubing to the patient as the last step.
- For disconnecting: shut the flowmeter (0 l/min position), remove the tubing from the patient, disconnect the humidifier if necessary, and then disconnect the flowmeter from the wall-mounted outlet or from the rapid outlet of the pressure regulator.

5. Symbols

	Do not use any grease, lubricant or oil.
	Refer to the instructions manual.
YYMMxxxx	Serial N° engraved on the device: YY: Year of manufacture, MM: Month of manufacture, xxxx: Unit number
	Manufacturer

6. Safety Instructions

- To prevent any risk of fire or explosion:
 - ✓ Do not smoke near the device.
 - ✓ Do not use near sources of ignition.
 - ✓ Do not allow the device to come into contact with greasy products.
 - ✓ Do not handle the device with oily hands.
 - ✓ Do not disconnect a device if it is under pressure.
- Do not restrict or block the outlet of the medical device.
- Do not use the flow outlet to operate any medical equipment.
- Influence of pressure: Important: check the supply pressure of the device. If the inlet pressure is not within the range indicated on the device, the flow may exceed the authorised standard range.
- Influence of temperature: Important: changes in ambient temperature (between 0°C and 40°C) affect the accuracy of the supplied flow rates.
- Do not disconnect the inlet or outlet fittings of the device.
- There is no indication showing the passage of gas on this type of device.
- The flow rate should be regulated on a snap-on position. Do not position the regulation knob between two flow rates.

7. Cleaning and disinfecting

- Disconnect the device from the gas supply. Clean the outside of the device with a clean, slightly damp cloth. Leave the device to dry. If disinfectants are used, check their compatibility with the plastic used (refer to the technical sheet).
- Do not immerse.
- Do not use a surface disinfectant.

8. Storage

Between -20 and +60°C in a clean and dry area.
If possible, store it in its packaging.

9. Maintenance

No preventive maintenance is required.

Only a person authorised by Technologie Médicale may dismantle the device.

Technologie Médicale recommends an annual check of the device:

- ✓ Verify that the device is properly connected to the wall-mounted outlet or to the outlet of the pressure regulator.
- ✓ Check that the gas flows correctly at every flow selection position.
- ✓ Ensure that there is no leakage when the knob is in the 0 position.

10. Warranty

The device is guaranteed for 1 year, parts and labour included. The warranty does not cover damage or accidents caused due to negligence, defective use, and lack of surveillance or maintenance.

Contractual service life: 10 years

Date of affixing the first **CE** marking: 2015

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تاريخ إلصاق العلامة © الأولى: 2015