C€0459

REGSON TM2

REGFLOW TM

DETREG TM

FR Notice d'instructions
EN Instructions for use
DE Bedienungsanleitung
ES Folleto de instrucciones
PT Manual de Instruções
IT Avvertenze e istruzioni
SV Bruksanvisning
NL Gebruiksaanwijzing
AR إرشادات الاستخدام

BG Инструкции за експлоатация
CS Návod k použití
EL Οδηγίες χρήσης
FI Käyttöohjeet
HU Használati útmutató
LT Naudojimo instrukcija
MT Ara n-nota ta' struzzjonijiet
PL Instrukcja obsługi
RO Manual de Instrucțiuni
SK Návod na použitie
SL Navodila za uporabo

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1. Presentation

Medical device meant to be connected on high pressure (HP) medical gas cylinders, to decrease the outlet pressure and supply a stable low pressure for use, and to adjust the outlet flow. The gauge indicates the pressure level of the gas cylinder.

User profile:

Users are medical staff: doctors, nurses.

For the REGSON TM2 and REGFLOW TM, users are also the home care technicians and the patient himself.

2. Versions

- Inlet connection and gas available: pin-index fitting in O₂ - Nut and nipple in : AFNOR in O₂, AIR and N₂O (for DETREG only) / DIN in O₂ and AIR / BS in O₂ / US in O₂ and AIR / UNI in O₂ / NORDIC in O₂ and AIR.
- REGSON TM2 model: pressure regulator delivered with a DEBSON TM2 pre-adjusted flowrates.
Available flowrates ranges: 0-1 l/min, 0-5 l/min, 0-15 l/min, 0-25 l/min, 0-50 l/min.
- REGFLOW model: pressure regulator delivered with a RTM3 flow meter with floating ball.
Available flowrates ranges: 0-1.5 l/min, 0-5 l/min, 0-15 l/min, 0-30 l/min.
- DETREG model: pressure regulator with AFNOR, BS or DIN quick-release connector and without flow meter.
- Threaded outlet for REGSON TM2 and REGFLOW: 12x125 M - 9/16" M - 1/4G M - 1/2" BS F.
- Quick-release connector in option for REGSON TM2 et REGFLOW TM.
- Injector in option.

All references and descriptions are available from the manufacturer.

3. Technical data

- In accordance with EN ISO 10524-1
- Unit: bar (pressure) and litre per minute (l/min for flow rate) for REGSON TM2 and REFFLOW TM.
- Pressure regulator designed with single stage piston.
- Built-in safety valve.
- Supply pressure: Max. 200 bar.
- Outlet pressure: 4.5 (± 0.5 bar).
- Working temperature: -20°C +60°C (-4°F +140°F).
- High pressure gauge accuracy: class 1.6.
- Flow rate accuracy:
REGFLOW TM: ±10% or ±0.5 l/min of the value, highest whichever.
REGSON TM2: ± 30 % ≤ 1,5 l/min ou ± 20 % > 1,5 l/min.
- On REGSON TM2, available flowrate between two positions.
- Serial number engraved on the body.

4. Instructions for use

Mounting the regulator on the cylinder:

- Make sure the device is undamaged, especially at the inlet connection.
- Slightly open, for a short time, the cylinder knob BEFORE connecting the pressure regulator, in order to remove impurities from the gas outlet.
Then turn it off.
- Check the presence and the good condition of the pressure regulator's input seal (nut and pin-index fitting model).
- Tightly screw, by hand, the inlet pin-index in the cylinder stopcock or the bull-nose on the cylinder.
- For the REGSON TM2, check that the position of the knob is on 0 l/min. For the REGFLOW, check that the flowmeter valve is closed.
- Gently open the cylinder stopcock.

DO NOT USE ANY WRENCH OR PLIERS

Flowmeter

- Screw, on the threaded outlet, a humidifier if necessary, or a nipple in order to connect tubing needed for therapy.
- Turn the adjustment knob or wheel on the front to the left until the desired pressure is displayed.
- Install the device to the patient (oxygen cannula for example) in last step.
- Equipped with a DEBSON TM2 barrel flow metre, the REGSON TM2 pressure regulator can be placed in any position.
- Equipped with a RTM3 ball flow metre, the REGFLOW pressure regulator must be placed in a vertical position.

Disassembling the regulator of the cylinder:

- First turn of the cylinder stopcock.
- The pressure regulator must be purged. The manometer must indicate 0 bar.
- For REGSON TM2 and REGFLOW TM, close the flowmeter.
- Remove the device from the patient and disassembly the humidifier if necessary.
- Unscrew by hand the pressure regulator from the cylinder.

Quick-release connector

Put the pressure regulator under pressure. Then, connect the standardised probe to the quick-release connector by checking the gas and connection type.

WHEN UNUSED, DO NOT LEAVE THE DEVICE UNDER PRESSURE

5. Symbols

	Use no oil		See instructions for use
AAMMxxxx	Serial number engraved on the device : AA : Year of manufacturing, MM : Month of manufacturing, xxxx : unique number		Manufacturer

6. Safety advice

IMPORTANT SAFETY REGULATIONS WHEN USING MEDICINAL OXYGEN CYLINDERS

EN

You are reminded that the operating personnel must be trained in the handling of gases:

- Never carry out multiple pressurisations successively.
- Never open the flow metre before opening the HP valve (it should always be set to 0 L/min beforehand).
- Do not close the HP valve with excessive torque (do not force).
- Do not use the threaded outlet to connect a pneumatic device.
- Do not stop up the outlets of the medical device.
- Influence of the temperature: Beware, the variation of the temperature has a real impact on the flow-rates accuracy.
- Do not disassembly the outlet and the inlet of the medical device.
- For the REGSON TM2: the knob should have been put in a real position. Do not try to adjust the knob between two positions of flow-rate. On this type of regulator, no indication shows the passage of gas.

Usual precautions with oxygen (O₂):

- **No smoking.**
- **Do not bring into contact with a flame.**
- **Do not lubricate.**
Do not coat the patients' face with greasy substances (Vaseline, ointments).
- Handle the equipment with clean hands, free from grease; the wearing of clean examination gloves (latex, nitrile rubber) is recommended.
- When opening, never stand in front of the cylinder valve output but always on the side opposite to the pressure regulator, behind the cylinder and taking a step back.
- Never expose the patient to the gas flow.
- Do not use aerosols (hairspray, deodorant etc.), solvents (alcohol, petrol etc.) on the equipment or in the vicinity.
- Check for leaks; in the event of a leak, close the valve. Never use a bottle presenting a sealing defect.
- Open the valve gradually.
- Never force the valve when opening.

You are reminded that the safe use of oxygen requires all the information on the label and in the instruction manual to be read carefully and any incident must be declared to the competent authorities.

7. Cleaning and disinfection

- Disconnect the regulator from the cylinder (see disassembly procedure).
- Use a disinfectant cleaner without alcohol for medical devices. Allow to dry before using again.
- Do not use floor cleaning products.
- Do not immerse

8. Storage

- Between -20 and +60°C (betw. 14 and 104°F) in a dry and clean place.
- Keep the device in its plastic bag as long as possible.

9. Maintenance

TECHNOLOGIE MEDICALE is the only allowed to operate on the device

The use of spare parts which are not in accordance with the manufacturer instructions exempts him from liability

Before each use

- Make sure the device is undamaged, due to a fall or shock, particularly on the gauge and the inlet connector
- Make sure there is a seal and its good condition on the inlet connector

At least once every 5 years

- Device overhaul by the manufacturer or a register provider: pressure regulator and flowmeter
- The overhaul deadline is displayed on the device

Technologie Médicale advice to control the device once a year:

- ✓ Check the good connection between the inlet and the cylinder valve.
- ✓ Check there is a flow on each position of the control knob (only for REGSON TM2).
- ✓ Check there is no leak when the control knob (for REGSON TM2) is at 0 l/min position or when the valve flowmeter is closed (for REGFLOW TM).

In case of fall, shock or leakage, have your device checked.

Spare parts	Use exclusively spare parts from TECHNOLOGIE MEDICALE
- Pin-index seals	
- Bull-nose seals	

Removal of the device does not have any special or unusual risks.

10. Guarantee

- 1 year warranty (parts and service), except in case of damage or accidents due to carelessness, misuse, bad supervision or maintenance.
- The closing of the regulation knob is guaranteed during 5 years without any leakage (REGFLOW TM).
- Contractual lifetime: 10 years.
- Date when **CE** marking first affixed: 2004

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- Datum prvního označení © 2004