
RTM3

CE 0459

FR	Notice d'utilisation
EN	Instruction for use
ES	Instrucciones de uso
DE	Bedienungsanweisung
IT	Istruzioni per l'uso
PT	Instruções de utilização
SV	Bruksanvisning
NL	Gebruikshandleiding
AR	تعليمات الاستخدام
BG	Инструкции за употреба
EL	Οδηγίες χρήσης
HU	Használati utasítás
LT	Naudojimosi instrukcija
MT	Struzzjonijiet għall-użu
PL	Instrukcja obsługi
RO	Instrucțiuni de utilizare
SK	Návod na používanie
SL	Navodila za uporabo
CS	Návod k použití
DA	Betjeningsvejledning
FI	Käyttöohjeet
LV	Lietošanas instrukcija
NO	Bruksanvisning
TR	Talimat bilgisi
RU	Инструкция по применению

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1. Presentation

The RTM3 and RTM3 DUO are medical devices intended to regulate the flow rate of a gas. They are designed to be connected to a pressurised wall-mounted gas supply or to the pressure outlet on a pressure regulator.

User profile:

Users are medical staff: doctors, nurses.

2. Versions

The RTM3 and RTM3 DUO are available in the following versions:

- Gas: Oxygen (O₂) - Medical air - Carbon dioxide (CO₂) and Carbogen (95% O₂ + 5% CO₂) on request.
- Pressure at device inlet: O₂ 4.5 bar (+/- 0.5 bar) Air 4.5 bar (+/- 0.5 bar)
- Probe: complies with the following standards: NF S 90-116 - BS 5682 - DIN 13260 - DISS - US Ohmeda Diamond - NORDIC - UNI....
- Direct connection or via a rail mounting (tube + crimped tube connector + rail clamp).
- Threaded outlet: 12x125 M - 9/16" M - 1/4G M - 1/2"BS F.
- Single or double outlet models.
- Flow ranges available:
 - ✓ O₂: 0-1.5 l/min; 0-5 l/min; 0-15 l/min; 0-30 l/min; 0-70 l/min.
 - ✓ Medical air: 0-5 l/min; 0-15 l/min; 0-30 l/min; 0-70 l/min.
 - ✓ CO₂: 0-12 l/min.
 - ✓ Carbogen: 0-15 l/min.

3. Technical data

- Complies with standard EN ISO 15002.
- Display units: litre per minute (l/min).
- Measurement precision: $\pm 10\%$ of the value indicated or ± 5 l/min, whichever is the highest.
- Maximum off-scale flow rate for ≤ 30 l/min version: 60 l/min $\pm 10\%$.
- Maximum off-scale flow rate for ≥ 30 l/min version: 90 l/min $\pm 10\%$.
- Inlet filters for ≤ 30 l/min version: mesh size of 35 μ m.
- Inlet filters for ≥ 30 l/min version: mesh size of 100 μ m.
- Serial number: with year, month of manufacture and individual number.
- The flowmeters provide (back-)pressure compensation (with the exception of the 0-1.5 l/min flowmeter).
- Ambient temperature for use: 23°C.
- RTM3 DUO \rightarrow 3 positions: gas flow through tubing nipple, gas flow through humidifier or OFF position.
- RTM3 DUO not available in 0-1.5 l/min version.

4. Instructions for use

- Check that the device is in good condition.
- Make sure that the regulating knob is closed (zero flow).
- Connect the device to the wall-mounted gas supply or to the pressure outlet on a pressure regulator, with the device positioned vertically.
- For the RTM3: If a humidifier is required, screw one onto the base of the device. Alternatively, screw on a tubing nipple to provide a connection for gas therapy equipment.
- For the RTM3 DUO: Screw a humidifier to the base of the device or connect gas therapy equipment directly to the tubing nipple. Select gas flow through the humidifier or through the tubing nipple by rotating the base of the device.
- Turn the regulating knob on the front of the device (anti-clockwise to increase the flow) until the required flow rate is displayed on the scale (**the value delivered is indicated by the middle of the ball**).
- For the RTM3 DUO: **Caution**, if the regulating knob is not closed, gas may escape at a high flow rate as you rotate the base to select an option other than the OFF position.
- Run the tubing up to the patient as the final step.
- When disconnecting: close the regulating knob, remove the tubing from the patient, disconnect the humidifier if one is fitted, then disconnect the flowmeter from the wall-mounted supply or from the pressure outlet on the pressure regulator.

5. Symbols

	Do not use grease, lubricant or oil
	Refer to the instructions
YYMMxxxx	Serial No. engraved on the device: YY: Year of manufacture, MM: Month of manufacture, xxxx: Individual number
	Manufacturer

6. Safety advice

- To prevent any risk of fire or explosion:
 - ✓ Do not smoke near the device.
 - ✓ Do not operate the device near a source of ignition.
 - ✓ Do not let the device come into contact with greasy or oily products.
 - ✓ Do not handle the device with greasy or oily hands.
 - ✓ Do not disassemble if the device is pressurised.
- Do not create an obstruction to restrict the outlet from the medical device.
- Effect of pressure: Caution: check the inlet pressure fed to the device. If this inlet pressure is not within the range indicated on the device, there is a danger that the flow rate delivered might be outside the permitted normative range.
- Effect of temperature: Caution: the precision of the flow rates delivered by the device varies as a function of the ambient temperature (between 0°C and 40°C).
- Do not disassemble the device's inlet or outlet connectors.
- Ensure that the device is always in a vertical position.
- Caution! Do not allow water from a connected humidifier to enter the flowmeter when it is removed from its normal operating position (for transport, cleaning etc.); always unscrew the humidifier before handling the device in any way.

7. Cleaning and disinfection

- Disconnect the device from the gas supply.
- Clean the outside of the device with a clean cloth, dampened with soapy water. Rinse the device then leave it to dry.
- Do not immerse in any liquid.
- Do not use an alcohol-based disinfection product.
- When using a disinfection product:
 - ✓ check that the product is compatible with the type of plastic used to make the device (polyamide, polypropylene and ABS) before using it.
 - ✓ check and comply with the product's directions for use (concentration, dilution, contact time, etc.)
- If necessary, contact the manufacturer if you have any question about these cleaning and disinfection instructions.

8. Storage

Store at between -20 and +60°C in a clean and dry place.

Store in its original packaging, if possible.

9. Maintenance

No preventive maintenance is required on the RTM3 during its service life.

Only individuals authorised by Technologie Médicale may disassemble the device.

Technologie Médicale recommends that the device should be inspected on a yearly basis:

- ✓ Check that the device can be connected correctly to a wall-mounted gas supply or to the pressure outlet on a pressure regulator.
- ✓ Check that the gas flows correctly at a range of different flow rates.
- ✓ Check that there is no leakage when the regulating knob is closed (i.e. zero flow).

10. Guarantee

The device is covered by a 1-year warranty (parts and labour), which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.

Contractual service life: 10 years.

Date when  marking first affixed: 2007

*In accordance with European Directive 93/42/EEC of June 14th 1993,
TECHNOLOGIE MEDICALE must be notified immediately of any incident or near-incident.*

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