C €0459

SOUPAPE DE JEANNERET

FR Notice d'instructions
EN Instructions for use
SK Návod na používanie

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1. Presentation

The SOUPAPE DE JEANNERET is a medical device for medical low suction.

User profile:

Users are medical staff: doctors, nurses.

2. Versions

The SOUPAPE DE JEANNERET is available with the following options:

- Connector: compliant with NF S 90-116 - BS 5682 - DIN 13260 - SS 87 524 30 - OHMEDA.
- Direct connection or via a rail mounting (tube + crimped tube holder + rail clamp).
- Outlet: 9mm diameter integrated nipple.
- Simple or double version.
- High of the water column: 60 or 90 cmH₂O.


3. Technical data

- Complies with ISO 10079-3 : 2014.
- Individual serial number: engraved on the device's body.
- 1 cm H₂O = 1 mbar = 0,1 KPa.
- Ambient temperature for use: 10 to 40°C.
- Atmospheric air pressure for storage and use: 800 – 1060 hPa.
- Relative humidity levels for storage and use: 0 to 100%.
- Minimum lighting levels for reading: 215 lx.

4. Instructions for use

Check the perfect tightness of the vacuum pipeline system.

- Check the valve **B** is well closed.
- Connect the device to the wall terminal unit.
- Fill the tube with water, up to a maximum level equal to 10 cm from the body: **X ≥ 10 cm**
- Unscrew the yellow nut **A** to release the graduated tube; adjust the requested suction level by setting the graduation in **cm** at the same level as water.
- Clamp the graduated tube by rescrewing the nut **A**.
- For the simple version, the graduated tube immersed in an **H** water level determines the suction level.
- For the double version, the water level inside the graduated tube at free flow, determines the suction level.
- Close the patient circuit and open the valve **B**, in order to obtain a bubbling as light and regular as possible (please note that a strong bubbling does not cause any variation of the water level **H**).
- Connect the device's outlet (nipple) to the collection jar using a suction tube with a minimum diameter of 6.3 mm.
- Close the valve **B** to stop the suction
- Proceed to a test before using on a patient.

5. Symbols

	Manufacturer
YYMMxxxx	Serial number engraved on the device : YY : Year of manufacturing, MM : Month of manufacturing, xxxx : serial number

6. Safety advice

- The SOUPAPE DE JEANNERET should only be used by qualified medical personnel.
- Always test the operation of the device prior to use.
- Make sure that the device is always fitted with an antibacterial filter in the suction circuit.
- Caution: the length of the tubing may have an impact on suction performance.
- Do not dispose of the antibacterial filters in a bin for domestic waste.
- Contact an authorised organisation to arrange for end-of-life disposal.
- The SOUPAPE DE JEANNERET is not compatible with use in an MRI environment.

7. Cleaning and disinfection

Always check compatibility of POLYCARBONATE with hospital disinfectants. Disinfectants containing raw FORMADEHYDE, GLYOXAL, AMPHOTERIC, GLUTARALDEHYDE or SUCCINIC ACID **are not compatible.**

- Dismantle components.
- Immerse the device in a pre-disinfecting solution, following manufacturer instructions and respecting the recommended dilution percentages and soaking times.
- Clean.
- Rinse abundantly and let the device in clean water during 1 minute, and then dry
- Sterilize or disinfect the device, following manufacturer instructions and respecting recommended dilution percentages and soaking times
- Put components back together

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Autoclavable **up to 121°C** or use cold sterilization.

Caution :

Polycarbonate is humidity and water-resistant. However, please note that a long immersion at a temperature over 60°C can lead to a decrease of the shock-resistance.

During the autoclaving, do not pile anything on the polycarbonate tubes.

8. Transport and storage

- Store between -20 and +60°C in a clean and dry place.
- Keep the instructions for use with the packaging for as long as possible.

9. Maintenance

Check every 1 to 3 years depending on level of use: Check the suction performance.

- When the device no longer connects up properly to the wall outlet, replace the connector. To replace the connector, use a suitable open-ended spanner. Check that the seal is in place before fitting the new connector.
- Caution: if the device has to undergo remedial maintenance work, an appropriate inspection and test must then be performed to check that the device can still be used safely. Record every maintenance operation.

In the event of damage or malfunction, return the device to the manufacturer or to an approved repairer.

Spare parts	Use exclusively
See the technical data sheet	TECHNOLOGIE MEDICALE spare parts

10. Guarantee

- The device is covered by a 1-year warranty (parts and labour), which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.
- Contractual service life: 10 years.
- Date of the first **CE** mark : 2000

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