VENTURI TM2

CE0459

FR	Notice d'instructions
EN	Instructions for use
DE	Bedienungsanleitung
ES	Folleto de instrucciones
PT	Manual de Instruções
IT	Avvertenze e istruzioni
BG	Инструкции за експлоатация
EL	Οδηγίες χρήσης
HU	Használati útmutató
PL	Instrukcja obsługi
RO	Manual de Instrucțiuni
LV	Lietošanas pamācība
LT	Naudojimo instrukcija
SR	Упутства за кориснике
SL	Navodila za uporabo
CS	Návod k použití
SV	Bruksanvisning
NL	Gebruiksaanwijzing
AR	إرشادات الاستخدام
DA	Brugsanvisning
SK	Návod na používanie

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1. Presentation
2. Versions
3. Technical data
4. Instructions for use
5. Symbols
6. Safety advice
7. Cleaning and disinfection
8. Transport and storage
9. Maintenance
10. Guarantee

1. Presentation

This medical device adjusts and measures the vacuum which it generates from a source of medical gas. It is used to aspirate fluids and mucus when a vacuum line is not available. It may be connected up to a source of oxygen or medical air. The VENTURI TM2 should be used in conjunction with a collection jar and a suction tube.

User profile:

Users are medical staff: doctors, nurses.

2. Versions

The VENTURI TM2 is available with the following options:

- Vacuum gauge display units available: 0-1000 mbar/hPa or 0-760 mmHg/0-100 kPa. 0-250 mbar/hPa or 0-200 mmHg/0-20 kPa
- Analogue or digital vacuum gauge.
- Connector: compliant with NF S 90-116 - BS 5682 - DIN 13260 - US Ohmeda Diamond - NORDIC - UNI - DISS.
- Direct connection or via a rail mounting (tube + crimped tube holder + rail clamp).
- Gases available: air or O₂.
- Outlet: integrated nipple, quick-fit safety jar, ½"G nipple or ½"G safety jar.

3. Technical data

- Complies with ISO 10079-3.
- Display units: millibar (mbar), hectoPascal (hPa) and millimetre of mercury (mmHg). 10 mbar = 10 hPa = 7.5 mmHg = 1 kPa.
- Measurement accuracy: ± 1.6 % of full scale (for analogue vacuum gauge) or ± 1 % of full scale (for digital vacuum gauge).
- Vacuum gauge can swivel through 90° (from - 45° to + 45°).
- Safety jar can swivel through 90° (from - 45° to + 45°).
- Antibacterial filter on front of device: 99.97% efficiency for 0.3 µm particles.
- Individual serial number: laser-engraved on the device's body.
- Safety valve fitted to protect device if exhaust becomes obstructed.
- Metal connection nipple, built into the body of the device.
- Operating pressure (supply): 3 to 6 bar.
- Maximum flow rate at 4.5 bar: 40 l/min with safety jar
- Maximum gas consumption at 4.5 bar: 60 l/min.
- Noise level with safety jar at maximum suction flow rate: 52 dB
- Ambient temperatures for use: 10 to 40°C.
- Atmospheric air pressure for storage and use: 800 – 1060 hPa.
- Relative humidity levels for storage and use: 0 to 100%.
- Minimum lighting levels for reading: 215 lx.

4. Instructions for use

VENTURI TM2

- Make sure that the device is undamaged and that the connector is compatible with the supply connection.
- Make sure that the regulation knob is in its closed position.
- Connect the device to the network's gas supply outlet (air or O₂).
- Connect the device's outlet (nipple or Ø 9 mm safety jar) to the collection jar using a suction tube with a minimum diameter of 6.3 mm.
- Turn the regulation knob on the front of the device anti-clockwise and set the vacuum level with the patient circuit closed off.
- To stop the suction, set the regulation knob to its closed position.

Quick-fit safety jar

- **Fitting:** insert and twist under the VENTURI TM2:
Line up the mark on the cover of the safety jar with the mark on the skirt under the body of the VENTURI TM2, then rotate the jar ninety degrees clockwise (moving the mark to the right).
- **Removing:** turn the safety jar to the left, line up the mark on the cover of the safety jar with that on the skirt under the body of the VENTURI TM2 then pull the jar downwards.

Plastic filter

- **Fitting:** push in fully until it clicks into place
- **Removing:** pull and twist at the same time.
- Dispose of the filter by observing the appropriate precautions.

5. Symbols

	Refer to the manual/instruction leaflet		Manufacturer
	Device is subject to separate collection		Vacuum adjustment indicator Vac = Vacuum

6. Safety advice

- The VENTURI TM2 should only be used by qualified medical personnel.
- Effect of supply gas pressure: caution, check the pressure delivered by the gas supply. This has an impact on the maximum vacuum level generated. The maximum vacuum levels and suction flow rates are only generated when the supply line pressure is 4.5 bar or greater.
- Always test the operation of the device prior to use.
- If the device is dropped, always make arrangements for the display accuracy to be checked.
- Never obstruct the device's air exhaust.
- Make sure that the needle (or digital display) indicates 0 when the device is not being used.
- Make sure that the device is always fitted with an antibacterial filter in the suction circuit.
- Caution: the length of the tubing may have an impact on suction performance.
- Caution: the measurement tolerance of the vacuum level display increases when the device is used outside its specified ranges of ambient temperature and pressure.
- Do not dispose of the antibacterial filters in a bin for domestic waste.
- Contact an authorised organisation to arrange for end-of-life disposal. The digital version of the VENTURI TM2 contains materials (electronic circuit and lithium batteries) that are toxic for the environment.
- When the source gas is oxygen (O₂), the use of lubricants that are not compatible with this gas may cause a fire or explosion.
- When the source gas is oxygen (O₂), ensure adequate ventilation of the room in which the VENTURI TM2 is used.
- Never dismantle the device when it is connected to a pressurised gas source.
- The VENTURI TM2 is not compatible with use in an MRI environment.

7. Cleaning and disinfection

Under normal conditions of use, the VENTURI TM2 does not require disinfection since the safety jar is protected by the filter on the front of the device. **Change the filter after each patient.** (See below for specific conditions).

VENTURI TM2

Use a disinfectant cleaner for medical devices. Leave to dry before using the device.

Caution:

- When using decontaminant products, make sure they are compatible with plastic (refer to the technical data sheet).
- Do not use surface decontaminants.
- Do not spray disinfectant cleaner directly onto the device. Use a cloth or a disposable wipe.
- Do not immerse the device in any product.

Safety jar

Since the safety jar is protected by an upstream filter on the front of the device, there is no need to sterilize it: simply replace the filter after each patient or, if the patient is admitted for a long period, replace as necessary after inspecting the filter to check the degree of fouling.

However, if the liquid accidentally overflows or if the filter perforates, then the safety jar must be disinfected or sterilized as follows:

1. Remove the safety jar and disassemble the cover from the jar.
2. Immerse safety jar parts in a pre-disinfecting solution (respecting the manufacturer's instructions concerning soaking, rinsing and drying times).
3. Clean, rinse thoroughly, soak for 1 minute in clean water and then dry.
4. Sterilize or disinfect in compliance with the instructions provided by the manufacturer of the particular product used.
5. Reassemble each part and refit the safety jar.

The safety jar can withstand 30 disinfection or sterilisation cycles.

The safety jar can be autoclaved at temperatures of up to 134°C.

8. Transport and storage

- To avoid generating a zero-offset error for the analogue vacuum gauge's needle, the device should be laid horizontal for transport.
- Store between -20 and +60°C in a clean and dry place.
- Keep the instructions for use with the packaging for as long as possible.

9. Maintenance

Check every 1 to 3 years depending on level of use: Check the accuracy of the vacuum gauge and the suction performance.

- Replace the silencer at least once a year or if liquid enters the device. Fouling of the silencer may result in a drop in suction performance. To remove the silencer, turn it anti-clockwise. To fit a new silencer, turn it clockwise, applying gently pressure.
- Replace the vacuum gauge if it displays inaccurate values. To replace the vacuum gauge, unscrew the two screws on the rear of the body.
- When the device no longer connects up properly to the wall outlet, replace the connector. To replace the connector, use a suitable open-ended spanner. Check that the seal is in place before fitting the new connector.
- Check that a seal is fitted to the outlet nipple. Its absence may cause a leak in the suction circuit and thus in a drop in performance.
- To replace the digital vacuum gauge's batteries, remove the vacuum gauge as described above. Remove the protective shells, then change the batteries using the foolproof connector.
- **Caution: Do not attempt to fit batteries using a connector other than those recommended by Technologie Medicale. The use of another connector could result in an explosion or fire.**
- Caution: if the device has to undergo remedial maintenance work, an appropriate inspection and test must then be performed to check that the device can still be used safely. Record every maintenance operation.

In the event of damage or malfunction, return the device to the manufacturer or to an approved repairer.

Spare parts	
- Silencer (Ref. 20195)	Use exclusively TECHNOLOGIE MEDICALE spare parts
- Complete safety jar, 100 ml (Ref. 18753)	
- Analogue vacuum gauge, 100 mbar, complete unit (Ref. 18739)	
- Analogue vacuum gauge, 760 mmHg complete unit (Ref. 19488)	
- Digital vacuum gauge, 100 mbar complete unit (Ref. 19776)	
- Digital vacuum gauge, 760 mmHg complete unit (Ref. 19779)	
- Lithium batteries for the digital vacuum gauge (Ref. 20094)	
- Tube of 10 filters (Ref. 11813)	
- Connector (contact us to determine the right Reference No.)	

10. Guarantee

- The device is covered by a 1-year warranty (parts and labour), which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.
- Contractual service life: 10 years.
- Date when **CE** marking first affixed: 2014

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