BLENDER TM2

OXYGENOTHERAPIE
MELANGEURS AIR-O₂
OXYGEN-THERAPY
AIR-O₂ BLENDERS

TECHNOLOGIE
MEDICALE

BLENDER TM2

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3 versions disponibles : Low Flow, High/Low Flow, High Flow.


The BLENDER TM2 is a medical device used to provide a mix of Oxygen and medical AIR to the patient (FiO_2).

It is intended to adjust the concentration of Oxygen (from 21% to 100%) in medical AIR.

It is connected to sources of pressurised gas on the wall or to pressure regulator's quick-release connectors using two hosepipes medical AIR and O_2 .

3 available versions: Low Flow, High/Low Flow, High Flow.

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Adjustment of the O_2 concentration (by %) (FiO_2). From 21% to 100 %. Percentage of O_2 precision : +/- 5 % of the full scale. In compliance with EN ISO 11195.

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
The RTM3-Leak features a unique and patented closed-circuit technology.

It guarantees the precision of the FiO_2 from 0 l/min on the Low Flow version and from 2 l/min on the High/Low Flow version. The leak is only activated when the flowmeter is open which allows significant gas savings (180 l/hour) in normal use and without any particular manipulation.

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The BLENDER TM2 mixer is fitted with protective ears for the adjustment knob.

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BLENDER TM2 is equipped with 1 outlet to connect the blender to an Oxygen analyzer.

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Audible alarm (80db to 1 m) in case of different pressure between Oxygen and Medical AIR.

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2 NIST standard inlet connectors to supply O_2 and medical AIR to the blender.

To be used with low pressure O_2 and medical AIR hosepipes with crimped probes: one NIST standard probe to connect to the blender + one standardized crimped probe (AFNOR, DIN, BS, US OHMEDA, DISS, UNI, NORDIC...) to be connected to a source of pressurized gas.

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The Air-O₂ Low Flow BLENDER TM2 is dedicated to Neonatal Intensive Care Unit (NICU).

- On the right the BLENDER TM2 is equipped with an integrated 5 l/min or 15 l/min or 30 l/min RTM3-Leak with a 3 l/min leak to guarantee the accuracy of the FiO₂ from 0 l/min.
- On the left, the BLENDER TM2 is equipped with an integrated 15 l/min RTM3-Leak or 30 l/min or simply an output to the DISS standard.

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The AIR-O₂ High/Low Flow BLENDER TM2 is used for high and low flow application. It is equipped with an integrated 5 l/min, 15 l/min or 30 l/min RTM3 Leak with a 10-12 l/min leak to guarantee the precision of the FiO₂ from 2 l/min and an integrated 15-70 l/min high flow RTM3.

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The AIR-O₂ High Flow BLENDER TM2 is mainly used in combination with a ventilator. It is equipped with a single integrated 15-70 l/min high flow RTM3.

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RTM3 flowmeter with floating ball O₂ 15 l/min mounted with DISS direct probe and outlet tubing nipple.

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Low pressure O₂ and medical AIR hosepipes with crimped probes: one NIST standard probe to connect to the blender + one standardized crimped probe (AFNOR, DIN, UNI, NORDIC...) to be connected to a source of pressurized gas.

Length: 3 meters as standard length. Other lengths available upon request (maximum 7 meters).

Lifetime: 10 years.

ACCESSOIRES DE FIXATION MOUNTING ACCESSORIES

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Complete trolley with BLENDER TM2 and O₂ analyzer.

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BLENDER TM2 and O₂ analyzer with adjustable alarm from 15% to 99% with non-permanent mount (thanks to an adapter).

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BLENDER TM2 and O₂ analyzer with adjustable alarm from 15% to 99% with permanent mount.

Aussi disponible / Also available

Distribué par / Distributed by

OXYGENOTHERAPIE
MELANGEURS AIR-O₂
OXYGEN-THERAPY
AIR-O₂ BLENDERS

TECHNOLOGIE
MEDICALE

101, rue Vaillant Couturier - 93130 NOISY-LE-SEC - FRANCE
Tél. : 33 (0)1 48 45 58 95
info@technologiemedicale.com - www.technologiemedicale.com