



SUCTION KITS & TROLLEYS

SUCTION
SUCTION KITS & TROLLEYS



**TECHNOLOGIE
MÉDICALE**



SUCTION KITS & TROLLEYS

The suction kits and trolleys enable to group on a mobile holder (either to carry or to push) the various products making up a suction set (suction devices, collection jars, catheter holders...).

SUCTION KITS:

Chrome-plated portable holder to put on the floor or to hang onto a rail thanks to a TM universal rail clamp in polycarbonate.

Dimensions: Height 370 mm x Width 320 mm x Depth 140 mm.

Ref. 10326: Portable suction kit only (without accessory).

Standard configuration:

- ▶ **One RVTM3** vacuum regulator 0-1000 mbar,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **One complete 2 L collection jar** in polysulfone with silicone ring lid,
- ▶ **One 2 L collection jar** in polysulfone (without lid),
- ▶ **One 3 m long silicone hose**, 7x13 mm, with one single-use vacuum stop.



CH01

Ref. 19428: RVTM3 0-1000 mbar with **AFNOR** probe.

Ref. 19429: RVTM3 0-1000 mbar with **BS** probe.

Ref. 19430: RVTM3 0-1000 mbar with **DIN** probe.

PORTABLE SUCTION TROLLEYS:

Rilsan-coated steel portable holder with 2 feet (at the front) and 2 wheels (at the back).

Dimensions: Height 750 mm x Width 400 mm x Depth 270 mm.

Ref. 11943: Portable suction trolley only (without accessory).

Standard configurations:

- ▶ **One RVTM3** vacuum regulator 0-1000 mbar,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **One complete 2 L collection jar** in polysulfone with silicone ring lid,
- ▶ **One 2 L collection jar** in polysulfone (without lid),
- ▶ **One 3 m long silicone hose**, 7x13 mm, with one single-use vacuum stop.



CH02

Ref. 19471: RVTM3 0-1000 mbar with **AFNOR** probe.

Ref. 19478: RVTM3 0-1000 mbar with **BS** probe.

Ref. 19479: RVTM3 0-1000 mbar with **DIN** probe.



- ▶ **One centralised vacuum inlet** with 2 vacuum supplies,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **Two complete 4 L collection jars** in polysulfone with pressure-fit lids and metal baskets,
- ▶ **Two 3 m long silicone hoses**, 8x14 mm, with two single-use vacuum stops.

CH03

Ref. 10242: with **AFNOR** probe.

*Other configurations available upon request.
Please contact us.*

RVTM3 also available in grey and/or with mmHg gauge on request.

MOBILE SUCTION TROLLEYS:

Chrome-plated mobile holder with a square-shaped column mounted on a 5 branch basis with 5 rotating wheels (including 2 wheels with brakes).
 Dimensions: Height 1140 mm x Width 500 mm x Depth 500 mm.
 Ref. 10215: Mobile suction trolley only (without accessory).

Standard configurations:

Single version

- ▶ One RVTM3 vacuum regulator 0-1000 mbar,
- ▶ One 3 m long connection to the central vacuum source with remote probe,
- ▶ One complete 2 L collection jar in polysulfone with silicone ring lid,
- ▶ One 2 L collection jar in polysulfone (without lid),
- ▶ One 3 m long silicone hose, 7x13 mm, with one single-use vacuum stop,
- ▶ One twin catheter holder, 40 cm long.

Twin version

- ▶ One twin RVTM3 vacuum regulator 1000-1000 mbar,
- ▶ One 3 m long connection to the central vacuum source with remote probe,
- ▶ Two complete 2 L collection jars in polysulfone with silicone ring lids,
- ▶ Two 3 m long silicone hoses, 7x13 mm, with two single-use vacuum stops,
- ▶ One twin catheter holder, 40 cm long.

CH04

Ref. 19434: Single version with **AFNOR** probe.



Ref. 19434	RVTM3 0-1000 mbar with AFNOR probe
Ref. 19443	RVTM3 0-600 mbar with AFNOR probe
Ref. 19455	RVTM3 0-1000 mbar with BS probe
Ref. 19446	Twin RVTM3 1000-1000 mbar with AFNOR probe
Ref. 19453	Twin RVTM3 600-600 mbar with AFNOR probe
Ref. 19449	Twin RVTM3 1000-600 mbar with AFNOR probe
Ref. 19457	Twin RVTM3 1000-1000 mbar with BS probe
Ref. 19462	Twin RVTM3 1000-1000 mbar with DIN probe
Ref. 11852	Switch: Accessory making it possible to switch from one collection jar to the other without stopping the suction process. To be used on a mobile suction trolley equipped with one RVTM3 vacuum regulator and two complete collection jars with lids.

- ▶ One centralised vacuum inlet with 2 vacuum supplies,
- ▶ One 3 m long connection to the central vacuum source with remote probe,
- ▶ Two complete 4 L collection jars in polysulfone with pressure-fit lids and metal baskets,
- ▶ Two 3 m long silicone hoses, 8x14 mm, with two single-use vacuum stops,
- ▶ One twin catheter holder, 40 cm long.



CH05

Ref. 10229: with **AFNOR** probe.
 Ref. 17115: with **DIN** probe.

CH06

It is possible to adapt various lengths of medical rail on our mobile suction trolley to fix different storage baskets (for catheters, gloves...).

Please contact us.



RVTM3 also available in grey and/or with mmHg gauge on request.

Other configurations available upon request.
 Please contact us.



Also available



Distributed by

SUCTION SUCTION KITS & TROLLEYS



TECHNOLOGIE MEDICALE

101, rue Vaillant Couturier B.P. 46 F-93136 NOISY-LE-SEC CEDEX France
Tel: 33 (0)1 48 45 58 95 - Fax: 33 (0)1 49 42 90 21 - 33 (0)1 48 45 29 00 - info@technologiemedicale.com
www.technologiemedicale.com