




SUCTION KITS & TROLLEYS



SUCTION
SUCTION KITS & TROLLEYS



**TECHNOLOGIE
MÉDICALE**



MOBILE SUCTION TROLLEY TUBALU

The mobile suction trolleys enable to group on a mobile holder the various devices making up a suction set (suction devices, collection jars, catheter holders...).



TUBALU trolley with central tube made of aluminium, and with 5 branch basis with 5 rotating wheels.

► **light and ultra-handly**

► **Customizable trolley**
to reply to all customer needs.

STEPS TO CUSTOMISE YOUR TROLLEY TUBALU

1

Choose your vacuum system

- ▶ RVTM3, VENTURI TM2 or JTM
- ▶ Single version or double version

Supplied in standard with 3 meters hose and a standardized crimped probe.



2

Choose your collection jars

- ▶ 1L, 2L or 4L jars
- ▶ Possibility to put up to 4 jars of 1L or 2L
- ▶ Single-use TM suction liners



3

Choose your accessories

- ▶ Single or double probe holder
- ▶ Vacuum stop
- ▶ ...



MOBILE SUCTION TROLLEYS CHROME-PLATED:

Chrome-plated mobile holder with a square-shaped column mounted on a 5 branch basis with 5 rotating wheels (including 2 wheels with brakes).

Dimensions: Height 1140 mm x Width 500 mm x Depth 500 mm.



CH01-E

- ▶ **Mobile suction trolley**
with RVTM3 0-1000 mbar with AFNOR probe.

- ▶ **One centralised vacuum inlet** with 2 vacuum supplies,
- ▶ **One 3 m long connection** to the central vacuum source with remote probe,
- ▶ **Two complete 4 L collection jars** in polysulfone with pressure-fit lids and metal baskets,
- ▶ **Two 3 m long silicone hoses**, 8x14 mm, with two single-use vacuum stops,
- ▶ **One twin catheter holder**, 40 cm long.



CH02-E

CH03-E

- ▶ **Contact us.**

It is possible to adapt various lengths of medical rail on our mobile suction trolley to different storage baskets (for catheters, gloves...).



RVTM3 also available in grey and/or with mmHg gauge on request.

Other configurations available upon request.
Please contact us.

SUCTION KITS & TROLLEYS

The suction kits and trolleys enable to group on a mobile holder (either to carry or to push) the various products making up a suction set (suction devices, collection jars, catheter holders...).

SUCTION KITS:

Chrome-plated portable holder to put on the floor or to hang onto a rail thanks to a TM universal rail clamp in polycarbonate.

Dimensions: Height 370 mm x Width 320 mm x Depth 140 mm.

Standard configurations:

- ▶ **One RVTM3** vacuum regulator 0-1000 mbar,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **One complete 2 L collection jar** in polysulfone with silicone ring lid,
- ▶ **One 2 L collection jar** in polysulfone (without lid),
- ▶ **One 3 m long silicone hose**, 7x13 mm, with one single-use vacuum stop.



CH04-E

PORTABLE SUCTION TROLLEYS:

Rilsan-coated steel portable holder with 2 feet (at the front) and 2 wheels (at the back).

Dimensions: Height 750 mm x Width 400 mm x Depth 270 mm.

Standard configurations:

- ▶ **One RVTM3** vacuum regulator 0-1000 mbar,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **One complete 2 L collection jar** in polysulfone with silicone ring lid,
- ▶ **One 2 L collection jar** in polysulfone (without lid),
- ▶ **One 3 m long silicone hose**, 7x13 mm, with one single-use vacuum stop.



CH05-E



- ▶ **One centralised vacuum inlet** with 2 vacuum supplies,
- ▶ **One 2 m long connection** to the central vacuum source with remote probe,
- ▶ **Two complete 4 L collection jars** in polysulfone with pressure-fit lids and metal baskets,
- ▶ **Two 3 m long silicone hoses**, 8x14 mm, with two single-use vacuum stops.

CH06-E

RVTM3 also available in grey and/or with mmHg gauge on request.

*Other configurations available upon request.
Please contact us.*



Also available



Distributed by

SUCTION SUCTION KITS & TROLLEYS



101, rue Vaillant Couturier - 93130 NOISY LE SEC FRANCE
Tel: +33 (0)1 48 45 58 95
info@technologiemedicale.com - www.technologiemedicale.com